

Product datasheet for **RR217309**

Grin1 (NM_001270605) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Grin1 (NM_001270605) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Grin1
Synonyms:	GluN1; NMDAR1; NR1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)


ORF Nucleotide Sequence:

>RR217309 representing NM_001270605
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGCACCATGCACCTGCTGACATTGCGCCTGCTTTTTCTGCTCCTTCGCCCGCCGCGCTGCGACC
 CCAAGATCGTCAACATCGGCGCGGTGCTGAGCAGCGCAAGCATGAACAGATGTTCCGCGAGGCAGTAAA
 CCAGGCCAATAAGCGACACGGCTCTTGGGAAGATACAGCTCAACGCCACTTCTGTACCCACAAGCCAAAC
 GCCATACAGATGGCCCTGTCAAGTGTGAGGACCTCATCTCTAGCCAGGTCTACGCTATCCTAGTTAGCC
 ACCCGCTACTCCAACGACCCTTCACTCCCACCCTGTCTCCTACACAGCTGGCTTCTACAGAATCCC
 TGTCTGGGACTGACTACCCGAATGTCCATCTACTCTGACAAGAGTATCCACCTGAGTTTCTTCGCACG
 GTGCCGCCCTACTCCACCAGTCCAGCGTCTGGTTTGAGATGATGCGAGTCTACAACCTGGAACCATCA
 TCCTGCTGGTCAAGGACGACCACGAGGGACGGCAGCGCAGAAGCGCTTGGAGACGTTGCTGGAGGAACG
 GGAGTCCAAGGCAGAGAAGGTGCTGCAGTTTGACCCAGGAACCAAGAATGTGACGGCTCTGCTGATGGAG
 GCCCGGGAACCTGGAGGCCCGGGTCAATCCTTTCTGCAAGCGAGGACGAGCTGCCACAGTGTACCGCG
 CAGCCGCAATGCTGAACATGACGGGCTCTGGGTACGTGTGGCTGGTGGGGAACGCGAGATCTCTGGGAA
 CGCCCTGCGCTACGCTCCTGATGGCATCATCGACTTCAGCTCATCAATGGCAAGAATGAGTCAGCCAC
 ATCAGTGACGCGGTGGCGTGGTGGCACAGGCAGTTCACGAACCTCTAGAGAAGGAGAATCACTGACC
 CACCGCGGGTTGCGTGGGCAACACCAACATCTGGAAGACAGGACCATTGTTCAAGAGGGTGTGATGTC
 TTCTAAGTATGCGGACGGAGTACTGGCCGTGTGGAATCAATGAGGATGGGACCGGAAGTTTGCCAAC
 TATAGTATCATGAACCTGCAGAACCGAAGTGGTCAAGTGGGCATCTACAATGGTACCCATGTCATCC
 CAAATGACAGGAAGATCATCTGGCCAGGAGGAGACAGAGAACTCGAGGATACCAGATGTCACCACG
 ACTAAAGATAGTGACAATCCACCAAGAGCCCTTGGTGTACGTCAAGCCCAATGAGTGTGGGACATGC
 AAAGAGGAGTTCACAGTCAATGGTGACCCAGTGAAGAAGGTGATCTGTACGGGGCCTAATGACACGTCCC
 CAGGCAGCCACGCCACACAGTGGCCAGTGTCTGCTATGGCTTCTGCATAGACCTGCTCATCAAGCTGGC
 GCGGACCATGAATTTTACCTATGAGGTGCACCTGGTGGCAGATGGCAAGTTTGGCACACAGGAGCGGGTA
 AACACAGCAACAAAAGGAGTGGAAACGGAATGATGGGCGAGCTACTCAGTGGCAAGCGGACATGATTG
 TGGCACCACTGACCATCAACAATGAGCGTGGCAGTACATAGAGTTCTCCAAGCCCTTCAAGTACCAGGG
 CCTGACCATTTTGGTCAAGAAGGAGATCCAGGAGCACACTGGACTCATTTATGCAGCCTTTTCAGAGC
 AACTGTGGTTGCTAGTAGGACTGTCAGTTCATGTGGTGGCTGTGATGCTGTACCTGCTGGACCGCTTCA
 GTCCTTTGGCCGATTCAAGGTGAACAGTGAGGAGGAGGAGGAAGATGCACTGACCCTGTCTCTGCCAT
 GTGGTTTTCTGGGGCGTCTGCTCAACTCCGGCATTGGGGAAGGTGCCCCCGGAGTTTCTCTGCACGT
 ATCCTAGGCATGGTGTGGGCTGGTTTCGCCATGATCATAGTGGCTTCTACACTGCCAAGTTGGCAGCTT
 TCCTGGTGTGGATCGGCCTGAGGAGCGCATCACGGGCATCAATGACCCAGGCTCAGAAACCCCTCAGA
 CAAGTTCATCTACGCAACTGTAAGCAGAGCTCCGTGGACATCTACTCCGGAGGCAGGTGGAGTTGAGT
 ACCATGTACCGGCACATGGAAAAACACAATTACGAGAGCGCAGCTGAGGCCATCCAGGCTGTGCGGGACA
 ACAAGCTGCACGCCTTTATCTGGGACTCGGCCGTGCTGGAGTTTGAGGCTTCACAGAAGTGGCAGTCTGGT
 GACCACGGGTGAGCTGTTCTCCGCTCAGGCTTTGGCATCGGCATGCGCAAGGACAGCCCTGGAAGCAG
 AACGTTTTCCCTGTCCATACTCAAGTCCCATGAGAATGGCTTCATGGAAGATCTGGATAAGACATGGGTTT
 GGTATCAGGAATGCGACTCCCGCAGCAATGCTCCTGCAACCCCTCACTTTTGAGAATGCGAGGGTCTT
 CATGCTGGTGGCTGGAGGCATCGTAGTGGGATTTTCTCATTTTCATTGAGATCGCCTACAAGCGACAC
 AAGGATGCCCGTAGGAAGCAGATGCAGCTGGCTTTTGCAGCCGTGAACGTGTGGAGGAAGAACCTGCAGA
 GCACCGGGGTGGACGCGGCGCTTTGCAAAACCAAAAAGACACAGTGTGCCGCGACGCGCTATTGAGAG
 GGAGGAGGCCAGCTGCAGCTGTGTTCCCGTCATAGGGAGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001270605

ORF Size: 2703 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001270605.1](#), [NP_001257534.1](#)

RefSeq Size: 4169 bp

RefSeq ORF: 2706 bp

Locus ID: 24408

UniProt ID: [P35439](#)

Cytogenetics: 3p13

MW: 101.3 kDa

Gene Summary: subunit of NMDA-preferring ionotropic glutamate receptors; may play a role in long term potentiation [RGD, Feb 2006]