

Product datasheet for SA013

Rat IgA (kappa) Rat Protein

Product data:

Product Type:	Native Proteins
Description:	Rat IgA kappa rat protein, 1 mg
Species:	Rat
Protein Source:	Myeloma
Concentration:	lot specific
Purity:	>95% Single band by Agarose Electrophoresis (5) and Immunoelectrophoresis (6).
Buffer:	State: Liquid purified monoclonal IgA Buffer System: PBS Preservative: 0.09% Sodium Azide
Preparation:	Liquid purified monoclonal IgA
Applications:	ELISA: 0.15-0.6 µg/ml.
Protein Description:	RAT IgA KAPPA, Clone IR22
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)