

## Product datasheet for SA015

### Rat IgG1 (lambda chain) Rat Protein

#### Product data:

Product Type:	Native Proteins
Description:	Rat IgG1 (lambda chain) rat protein, 1 mg
Species:	Rat
Protein Source:	Myeloma
Concentration:	lot specific
Purity:	>95% pure
Buffer:	State: Liquid IgG, purified by affinity chromatography Buffer System: Phosphate buffered saline, 0.09% Sodium Azide
Preparation:	Liquid IgG, purified by affinity chromatography
Applications:	ELISA.
Protein Description:	IgG1 Lambda Myeloma Protein, Clone IR31. The antigen specificity of this product is unknown.
Storage:	Store the antigen undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Centrifuge vial before opening. Avoid repeated freezing and thawing. Should this product contain a precipitate we recommend microcentrifugation before use.
Stability:	Shelf life: one year from despatch.
Summary:	SA015 is a spontaneously arising rat IgG1 myeloma.



[View online »](#)