

Product datasheet for SA021

Rat IgG2a Rat Protein

Product data:

Product Type:	Native Proteins
Description:	Rat IgG2a rat protein, 1 mg
Species:	Rat
Protein Source:	Myeloma
Concentration:	lot specific
Purity:	>95%
Buffer:	State: Liquid Buffer System: PBS containing 0.09% Sodium Azide
Preparation:	Liquid
Applications:	ELISA.
Protein Description:	Purified myeloma protein prepared by affinity chromatography.
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)