

Product datasheet for SA022

Rat IgE Kappa Rat Protein

Product data:

Product Type:	Native Proteins
Description:	Rat IgE Kappa rat protein, 0.1 mg
Species:	Rat
Protein Source:	Myeloma
Concentration:	lot specific
Purity:	>95% pure. Single band by Agarose Electrophoresis (Ref. 5) and Immunoelectrophoresis (Ref. 6) Preparation: Affinity Chromatography (Ref. 4)
Buffer:	State: Liquid purified IgE (myeloma IR162) Buffer System: PBS, pH 7.2 containing 0.09% Sodium Azide as preservative
Preparation:	Liquid purified IgE (myeloma IR162)
Applications:	ELISA: 0.1-1.25 µg/ml.
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)