

Product datasheet for SA023

Rat IgG2b (kappa chain) Rat Protein

Product data:

Product Type:	Native Proteins
Description:	Rat IgG2b (kappa chain) rat protein, 1 mg
Species:	Rat
Protein Source:	Myeloma
Concentration:	lot specific
Purity:	>95% Single band by agarose electrophoresis and immunoelectrophoresis
Buffer:	State: Liquid Buffer System: PBS, pH 7.2 containing 0.09% Sodium Azide
Preparation:	Liquid
Applications:	ELISA: 0.3 µg/ml.
Protein Description:	Purified myeloma IgG2b KAPPA (Clone:IR863)
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)