

## Product datasheet for SA025

### Rat IgM Rat Protein

#### Product data:

Product Type:	Native Proteins
Description:	Rat IgM rat protein, 1 mg
Species:	Rat
Protein Source:	Myeloma
Concentration:	lot specific
Purity:	>95%
Buffer:	State: Liquid Buffer System: PBS containing 0.09% Sodium Azide
Preparation:	Liquid
Applications:	ELISA: 0.1 - 1.0µg/ml.
Protein Description:	Purified IgM prepared by affinity chromatography (Clone: IR473).
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Rat Immunoglobulin M
Protein Families:	ELISA: 0.1 - 1.0µg/ml.



[View online »](#)