

Product datasheet for SA6013

SNAP25 (1-226, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	SNAP25 (1-226, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MAEDADMRNE LEEMQRRADQ LADESLESTR RMLQLVEESK DAGIRTLVML DEQGEQLERI EEGMDQINKD MKEAEKNLTD LGKFCGLCVC PCNKLKSSDA YKKAWGNNQD GVASQPARV VDREQMAIS GGFIRRVNTD ARENEMDENL EQVSGIIGNL RHMALDMGNE IDTQNRQIDR IMEKADSNKT RIDEANQRAT KMLGSG
Tag:	His-tag
Predicted MW:	25.4 kDa
Concentration:	lot specific
Purity:	>90% by SDS-PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 50 mM HEPES pH 7.4, 5 mM DTT, 100 mM NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human SNAP25, fused to His tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001309831
Locus ID:	6616
UniProt ID:	P60880
Cytogenetics:	20p12.2
Synonyms:	bA416N4.2; CMS18; dJ1068F16.2; RIC-4; RIC4; SEC9; SNAP; SNAP-25; SUP



[View online »](#)

Summary:

Synaptic vesicle membrane docking and fusion is mediated by SNAREs (soluble N-ethylmaleimide-sensitive factor attachment protein receptors) located on the vesicle membrane (v-SNAREs) and the target membrane (t-SNAREs). The assembled v-SNARE/t-SNARE complex consists of a bundle of four helices, one of which is supplied by v-SNARE and the other three by t-SNARE. For t-SNAREs on the plasma membrane, the protein syntaxin supplies one helix and the protein encoded by this gene contributes the other two. Therefore, this gene product is a presynaptic plasma membrane protein involved in the regulation of neurotransmitter release. Two alternative transcript variants encoding different protein isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome

Protein Pathways:

SNARE interactions in vesicular transport

Product images: