

Product datasheet for SA6031X

HSPB1 / HSP27 Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	HSPB1 / HSP27 human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MTERRVPFSL LRGPSWDPFR DWYPHSRLFD QAFGLPRLPE EWSQWLGSS WPGYVRPLPP AAIESPAVAA PAYSRALSQR LSSGVSEIRH TADRWRVSLD VNHFAPELDT VKTKDGVVEI TGKHEERQDE HGYISRCFTR KYTLPPGVDP TQVSSLSPE GTLTVEAPMP KLATQSNEIT IPVTFESRAQ LGGPEAAKSD ETAAK
Predicted MW:	22.7 kDa
Concentration:	lot specific
Purity:	>95% pure by SDS-PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Hepes (pH 7.5) containing 1 mM DTT, 100 mM KCl
Endotoxin:	< 1.0 EU per 1µg of protein (determined by LAL method)
Preparation:	Liquid purified protein
Protein Description:	Recombinant Human Hsp27 was overexpressed in <i>E. coli</i> and purified by conventional chromatography.
Storage:	Store the antigen at 2-8°C for 1-2 weeks or (in aliquots) at -20°C to -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001531
Locus ID:	3315
UniProt ID:	P04792 , V9HW43
Cytogenetics:	7q11.23
Synonyms:	CMT2F; HEL-S-102; HMN2B; HS.76067; Hsp25; HSP27; HSP28; SRP27



[View online »](#)

Summary:

This gene encodes a member of the small heat shock protein (HSP20) family of proteins. In response to environmental stress, the encoded protein translocates from the cytoplasm to the nucleus and functions as a molecular chaperone that promotes the correct folding of other proteins. This protein plays an important role in the differentiation of a wide variety of cell types. Expression of this gene is correlated with poor clinical outcome in multiple human cancers, and the encoded protein may promote cancer cell proliferation and metastasis, while protecting cancer cells from apoptosis. Mutations in this gene have been identified in human patients with Charcot-Marie-Tooth disease and distal hereditary motor neuropathy. [provided by RefSeq, Aug 2017]

Protein Pathways:

MAPK signaling pathway, VEGF signaling pathway

Product images: