

Product datasheet for SA6053

MIF (1-115, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	MIF (1-115, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MRGSHHHHHH</u> GMASMTGGQQ MGRDLYDDDD KDRWGSMPMF IVNTNVPRAS VPDGFLSELT QQLAQATGKP PQYIAVHVVP DQLMAFGGSS EPCALCSLHS IGKIGGAQNR SYSKLLCGLL AERLRISPDR VYINYYDMNA ANVGWNNSTF A
Tag:	His-tag
Predicted MW:	16.6 kDa
Concentration:	lot specific
Purity:	>95% by SDS-PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: PBS, pH 7.4
Endotoxin:	< 1.0 EU per 1 µg of protein (determined by LAL method)
Preparation:	Liquid purified protein
Protein Description:	Recombinant human MIF, fused to His-tag at N-terminus, was cloned into an E. coli expression vector and was purified to apparent homogeneity by using conventional column chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_002406</u>
Locus ID:	4282
UniProt ID:	<u>P14174</u> , <u>I4AY87</u>
Cytogenetics:	22q11.23
Synonyms:	GIF; GLIF; MMIF



[View online »](#)

Summary:

This gene encodes a lymphokine involved in cell-mediated immunity, immunoregulation, and inflammation. It plays a role in the regulation of macrophage function in host defense through the suppression of anti-inflammatory effects of glucocorticoids. This lymphokine and the JAB1 protein form a complex in the cytosol near the peripheral plasma membrane, which may indicate an additional role in integrin signaling pathways. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome

Protein Pathways:

Phenylalanine metabolism, Tyrosine metabolism

Product images: