

Product datasheet for SC101293

TMEM167A (NM_174909) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM167A (NM_174909) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM167A
Synonyms:	TMEM167
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF within SC101293 sequence for NM_174909 edited (data generated by NextGen Sequencing) ATGTCTGCCATTTTCAATTTTCAGAGTCTATTGACTGTAATCTTGCTGCTTATATGTACC TGTGCTTATATTTCGATCCTTGGCACCCAGCCTCCTGGACAGAAATAAACTGGATTGTTG GGTATATTTTGGAAAGTGTGCCAGAATTGGTGAACGGAAGAGTCCTTATGTTGCAGTATGC TGTATAGTAATGGCCTTCAGCATCCTCTTCATACAGTAG Clone variation with respect to NM_174909.4 >OriGene 5' read for NM_174909 unedited GCTGCGGCATCGCCGGCCTCGCTATGTCTGCCATTTTCAATTTTCAGAGTCTATTGACTG TAATCTTGCTGCTTATATGTACCTGTGCTTATATTTCGATCCTTGGCACCCAGCCTCCTGG ACAGAAATAAACTGGATTGTTGGGTATATTTTGGAAAGTGTGCCAGAATTGGTGAACGGA AGAGTCCTTATGTTGCAGTATGCTGTATAGTAATGGCCTTCAGCATCCTCTTCATACAGT AGCTGGGGAATGCCAGAATGTAGTTGCCATCAGATTTGATTGTGAACAAGGACTGACT GCAGAAATAATGGAAGGATGTTAACTCTTTATCTCCGAACATTGAATGAGATAAAT TTCCAGATGCTGTTCTCTATTTTAATGTTATTGGACCAATGTTCTGTATAACAATTAAG ATGTAACCATTTAATAGTCTGTAAACAATCAACCTCAGTACTGTCACTACAATATTACATT CTGCAAAATGTNATTCTGTTGTATCAGATACAAAATTTTAGTGAGGTATCTCTAANGCACA TAGTAGAAAACAAAATTGGTTAATTACTCAAGTTCCTTTCACTGNGATTGGAAATGATT TAATCTTTATAGAATGAAAACCTTTTTTNGGACTAACCTTTTTTATTAATGGCTCAAT TTTGNGTTGATAAGGATTGCATTAATATTTNAATAGGGGCTTGCTTTTCTCTGGGCC ACCATTTTGGATCATTAAACAAAACCTCTACTCNTAACNACTCTANTTTATGACAAG NNTTAAAAATTTTAAACAAGCTTATGCAG
5' Read Nucleotide Sequence:	
Restriction Sites:	NotI-NotI
ACCN:	NM_174909


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_174909.1</u> , <u>NP_777569.1</u>
RefSeq Size:	726 bp
RefSeq ORF:	219 bp
Locus ID:	153339
UniProt ID:	<u>Q8TBQ9</u>
Cytogenetics:	5q14.2
Protein Families:	Transmembrane
Gene Summary:	Involved in the early part of the secretory pathway.[UniProtKB/Swiss-Prot Function]