

## Product datasheet for **SC103754**

### ZNF573 (BC029468) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF573 (BC029468) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF573
Synonyms:	FLJ30921
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL4</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for BC029468, the custom clone sequence may differ by one or more nucleotides

```
ATGACCTGTGAATTAGTGACATTCAGGGATGTGGCCATAGACTTCTCTCGGCAGGAGTGGGAATACCTGG
ACCCTAATCAGAGGGACTTATACAGGGATGTGATGTTGGAGAACTATAGAAACCTGGTATCACTGGGAGG
ACATTCCATTTCTAAACCAGTTGTGGTTGATTTACTGGAGCGAGGAAAAGAGCCCTGGATGATTTTGAGG
GAAGAAACACAGTTCACAGATTTGGATTTACAGTGTGAGATAATCAGCTACATAGAAGTACCCACTATG
AAACAGATATATCCTCTACACAACCTCAGAGCATATATAAGAGAGAGAACTCTATGAATGA
```



[View online »](#)

<b>5' Read Nucleotide Sequence:</b>	>OriGene 5' read for BC029468 unedited CAACATTTTGTAAATACGACTCACTTATAGGGCGGCCGCGATTTCGGCAGGAGTTCCCTTC CGTTTTTCGCGTGGTTCTTTTGAAGCTCTGGATTCTCTGGAGTTTGAATGTTTCCAGTAT TGGAAACCCACCAAGTAGGACTGATCAGGTCTTACAATTCTAAAACCATGACCTGTTTTTC AGGAATTAGTGACATTCAGGGATGTGGCCATAGACTTCTCTCGGCAGGAGTGGGAATACC TGGACCCTAATCAGAGGGACTTATACAGGGATGTGATGTTGGAGAATATAGAAACCTGG TATCACTGGGAGGACATCCATTTCTAAACCAGTTGTGGTTGATTTACTGGAGCGAGGAA AAGAGCCCTGGATGATTTTGAGGGAAGAAACACAGTTCACAGATTGGATTACAGTGTG AGATAATCAGCTACATAGAAGTACCCACTTATGAAACAGATATATCCTCTACACAACCTC AGAGCATATATAAGAGAGAGAACTCTATGAATGTAAGAAATGTCAGAAGAAATTTAGTA GTGGTTATCAACTTATTCTACATCACAGGTTTCATGTCATTGAGAGACCCTATGAATGCA AAGAGTGTGGGAAGAAGCTTCGTAGTGGCTATCAACTACTCTACATCAAAGATTTTCATA CTGGTGAGAAACCCTATGAATGTACCGAATGTGGGGAAGAAGCTTTAGAAGTGGGTATCAG CTGACTGTGCATCAGAGATTTTCACTGGTGAGAAAACCTATGAATGTACGCAGTGTGGG AAGGCCTTTATATATGCCTCAAACCTTGTCAACATGAGAAGAATCACACTGGGTGGCAA CCCGCAGAATGTCCGGAGTGTGGGAAGGCCTCTTATTACAGGGGACATCTTATAAATCT TATAAAAGCCTACCTGCCGAAAACCATATAAAGTACAGAATGTGGGAAGACCTTTATAGG CGGCCAAAATGTGTGACGA
<b>Restriction Sites:</b>	NotI-NotI
<b>ACCN:</b>	BC029468
<b>Insert Size:</b>	2500 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">BC029468.1</a>
<b>RefSeq Size:</b>	720 bp
<b>RefSeq ORF:</b>	608 bp
<b>Locus ID:</b>	126231
<b>Cytogenetics:</b>	19q13.12
<b>Domains:</b>	KRAB
<b>Protein Families:</b>	Transcription Factors

**Gene Summary:** May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]