

Product datasheet for **SC104931**

COMMD7 (AK027893) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	COMMD7 (AK027893) Human Untagged Clone
Tag:	Tag Free
Symbol:	COMMD7
Synonyms:	C20orf92; dj1085F17.3
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Cell Selection:	None
Fully Sequenced ORF:	>NCBI ORF sequence for AK027893, the custom clone sequence may differ by one or more nucleotides

```
ATGGGCCGCTGCACTGCACTGAGGACCCGGTGCCGGAGGCCGTGGCGGCGACATGCAGCAGCTGAACC
AGCTGGGCGCGCAGCAGTTCTCAGCCCTGACAGAGGTGCTTTTCCACTTCCTAACTGAGCCAAAAGAGGT
GGAAAGATTTCTGGCTCAGCTCTCTGAATTTGCCACCACCAATCAGATCAGTCTTGGCTCCCTCAGAAGC
ATCGTGAAAAGCCTCCTTCTGGTTCCAAATGGTGCTTTGAAGAAGAGTCTCACAGCCAAGCAGGTCCAGG
CGGATTTCAACTCTGGGTCTTAGTGAGGAGAAAGCCACTTACTTTTCTGAAAAGTGGAAGCAGAATGC
TCCCACCCTTGCTCGATGGGCCATAGGTCAGACTCTGATGATTAACCAGCTCATAGATATGGAGTGGAAA
TTTGGAGTGACATCTGGGAGCAGCGAATTGGAGAAAGTGGGAAGTATATTTTACAACATAAAGTTGGTGG
TTAAGAAAGGAAATCAAACCGAAAATGTGTATATAGAATTAACCTTGCCTCAGTTCTACAGCTTCTGCA
CGAGATGGAGCGAGTCAGAACCAGCATGGAGTGTCTGCTGA
```



[View online »](#)

