

Product datasheet for **SC104969**

C16orf13 (METTL26) (NM_032366) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C16orf13 (METTL26) (NM_032366) Human Untagged Clone
Tag:	Tag Free
Symbol:	C16orf13
Synonyms:	C16orf13; JFP2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_032366, the custom clone sequence may differ by one or more nucleotides

```

ATGCTGGTGGCGCGGCCGCGGAGCGGAACAAGGATCCCATCTTGACGTGCTGCGGCAGTACCTGGATC
CGGCCAGCGTGGCGTCCGCGTCTCGAGGTGGCCTCGGGCTCCGGCCAGCACGCAGCGCACTTCGCGCG
GGCCTTCCCCCTGGCCGAGTGGCAGCCGTGGACGTGGACCAGCGCTGCCTGGACAGCATCGCGGCCACC
ACGCAAGCCCAGGGCCTGACCAACGTGAAGGCCCCGCTACACCTGGACGTGACGTGGGGCTGGGAGCACT
GGGGCGGGATCCTGCCACAGTCGCTGGACCTGTTGCTCTGCATCAACATGGCCCATGTGACCCCCCTGCG
CTGCACGGAGGGGCTCTCAGAGCAGCAGGACACCTGCTCAAACCCAGGGCCCTGCTCATCACCTACGGG
CCCTATGCCATCAATGGGAAGATCTCCCCCAGAGCAACGTGGACTTTGACCTGATGCTCAGATGCAGGA
ACCCAGAATGGGGGCTTCGGGACACAGCCCTCCTGGAGGACCTGGGAAAGGCCAGTGGCCTGCTCCTGGA
GAGGATGGTGGACATGCCAGCCAACAACAAATGCCTGATCTCCGAAAACTAA
  
```


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for NM_032366 unedited CACGAGGCCGCGTCCGCGCGGGGAGTTGTTGCTGCCGCGATGCTGGTGGCGGCGGCCGCGGAGCGGAACAAGGATCCCATCTTGACAGTGTGCGGCAGTACCTGGATCCGGCCAGCGTGGCGTCCGCGTCTCGAGGTGGCTCGGGCTCCGGCCAGCACGACGCGCACTTCGCGCGGGCCTTCCCCCTGGCCGAGTGGCAGCCGTCGGACGTGGACCAGCGCTGCCTGGACAGCATCGCGGCCACCAAGCAAGCCAGGGCCTGACCAACGTGAAGGCCCGCTACACCTGGACGTGACGTGGGGCTGGGAGCACTGGGGCGGGATCCTGCCACAGTCGCTGGACCTGTTGCTCTGCATCAACATGGCCCATGTACAGCCCTGCGCTGCACGGAGGGGCTCTTCAGAGCAGCAGGACACCTGCTCAAACCCAGGGCCCTGCTCATCACCTACGGGCCCTATGCCATCAATGGGAAGATCTCCCCCAGAGCAACGTGGACTTTGACCTGATGCTCAGATGCAGGAACCCAGAATGGGGCTTCGGGACACAGCCCTCTGGAGGACCTGGGAAAGGCCAGTGGCCTGCTCCTGGAGAGGATGGTGGACATGCCAGCCAACAACAAATGCCTGATCTTCCGGAAAACTAAGCCCCTNCTTACCCNCGCACACCTGCATCCCTGCCGGGAGCTCTGTGAGGCACGAACCTGCCTCCCT
3' Read Nucleotide Sequence:	>OriGene 3' read for NM_032366 unedited CTAGACAAACTNGTTTGGCAGCTGGGCCCTTCCGCCGCGGCTGAGACACAGACTGGTGGGGCTGTGCTCCACAAGGTCGGCCTACGGAGGCAGGGTTCTGCTCACAGAGCCTCCGGCAGGGATGCAGGTGTGCGGGGTGAAGGAGGGGCTTAGTTTTTCCGGAAGATCAGGCATTTGTTGTTGGCTGGCATGTCCACCATCCTCTTAGGAGCAGGCCACTGGCCTTCCAGGTCCTCCAGGAGGGTGTGTCCGAAGCCCCATTCTGGGTTCTGCATCTGAGCATCAGGTCAAAGTACACGTTGCTCTGGGGGAGATCTTCCATTGATGGCATAGGGCCCGTACGTGATGAGCAGGGCCCTGGGTTTGAGCAGGTGTACTGCTGCTCTGAAGAGCCCTCCGTGCAGCGCAGGGGGCTGACATGGCCCATGTTGATGCACAGCAACAGGTTACGCGATTGTGGCTGGATCCCGCCCCAGTGCTTTCATTCCACGCTACGTTACGGTGTATCCGGGACCTTACGCTTGTGAGGCCCTTGTCTTGTTCGGTGGTTCCCCCGCCGTTTCATGTCTTCTGGTCCCCCTCCC CGGTTGCTTTTTTTCTCCGGTCACGCCTCTCTCCCATAGTCTTCCCCTCGTGTCCCTTATCTCCTGACCTCTTTTTGTTGTTACTCTCCTCATTGTTTCGCCTCCCCTCCCATTCTTCTCCCTCCCCTCCCGTGTCTTCTTCCCCCCTCCTCCTCTCCCTCTTTTTTGCCAACCTTCCCCCTCCCCCTTTATTTCTTCATCCCCTCTGTCTTCCCGCCCGCGTCTCTCCACGTTTTCTCCTTTCTCCCCCTTTCCGCTTCTCCCCTCTATTCTCCCCGTCCCGT TTCTCTTCCCCCTCCCACCCACACCCATGTATCGCGTTGTTTCCGTCCGATTCTTCCTTCTCTTCTCCCTCGTCCCCGCTATCGACTCCGCCCTGTCTCCGTTGTTTGTCTCGATCTCCCACTACCCCTTCTCTATTCAT
Restriction Sites:	NotI-NotI
ACCN:	NM_032366
Insert Size:	730 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032366.1](#), [NP_115742.1](#)

RefSeq Size: 861 bp

RefSeq ORF: 861 bp

Locus ID: 84326

UniProt ID: [Q96S19](#)

Cytogenetics: 16p13.3