

Product datasheet for **SC105796**

LOC400627 (AF289562) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC400627 (AF289562) Human Untagged Clone
Tag:	Tag Free
Symbol:	LOC400627
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for AF289562, the custom clone sequence may differ by one or more nucleotides

ATGGAAGGGCTCTTGGGATGGGTCTTGGGCTGAGCCTTTGGCAAGGCCTCCACCCTGGAGCTGAGCCAG
CAGAGAGTGGGGCTGACTCCATGAGGCCAGACTCTCGTGGGTATTCTGGTGGCCTGGGCACCCTGTCTT
TCCTGAACCCCTGCCTGCCTAATCTAAGAACTGCTCCCCTTGTAGCTCTGTCTCCTCTGTTGGCCTAGT
TTCTTCATCCATGGGTGGGCACAGCAGAAGCTGTCTCACCTGCTCCATGGTGAACGGGATTTGCAAGGAT
TTGACTCTGGCCATCCAGGAAGATGGAGTCTTTGACCAGGCCTGGCCCTGA

5' Read Nucleotide Sequence:

>OriGene 5' read for AF289562 unedited
NGTTCACATTTTGTATACGACTCACTATAGGCGGCCGGAATCGGCACGAGGGTTTGGAG
GGGAGAGAGTGGTGCACAGGGACTGGAGGTGGGGGAGGGGTGGGGTTGCCATGGAAGGG
CTCTTGGGATGGGTCTTGGGCTGAGCCTTTGGCAAGGCCTCCACCCTGGAGCTGAGCCAG
CAGAGAGTGGGGCTGACTCCATGAGGCCAGACTCTCGTGGGTATTCTGGTGGCCTGGGC
ACCCTGTCTTTCCTGAACCCCTGCCTGCCTAATCTAAGAACTGCTCCCCTTGTAGCTCTGT
CTCTCCTCTGTTGGCCTAGTTTCTTCATCCATGGGTGGGCACAGCAGAAGCTGTCTCAC
TGCTCCATGGTGAACGGGATTTGCAAGGATTTGACTCTGGCCATCCAGGAAGATGGAGT
CTTTGACCAGGCCTGGCCCTGAAGCAGACTTCTCCTCTGAGCCTGCTGGGGAAGCAG
AATTCTCGATGCTGGGTGTCATGGTCTACCAGGCCCTCCCTCGCCAGGCCACTTGTGT
TGAGCCTGGAGCAGAACCCAGATCTGTAGTCCGCTGTGAGCTGGGGAGTCTCCGTT
TCCGTCTGTTTTCATCCATCAGCTGTGTGGTGTGAGCAAGTCACTTAACCTCTTGAGG
CCTCCATGTCCTCATCTGCAAAATGGAGACACACATATCTGCGTCCCGATGACCACCAG
TAGCTGGGATGATACGTCGACTGCAGAGTGACCACTTTGTGTAGGCGTCCCTGTGGATGG
TTGTCAGTCAAAGGAGTGAGTCANATGATGAGAGGGGAGGACAGCANAGCAGCAGNTGTT
TCTTCAGTCTCGCCGGGACTCCAGGATCGCCATGTTCTGCAGCCCCGCCTCTACTGTCT
CCCTGGTGGGACCTGAACCTCAACTGCTTAGTTTAAAGTCCG



[View online »](#)

Restriction Sites:	NotI-NotI
ACCN:	AF289562
Insert Size:	1600 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	AF289562.1 , AAL55746.1
RefSeq Size:	2131 bp
RefSeq ORF:	333 bp
Locus ID:	400627
Cytogenetics:	17q25.3