

Product datasheet for **SC105842**

GIDRP88 (R3HCC1L) (AF050198) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GIDRP88 (R3HCC1L) (AF050198) Human Untagged Clone
Tag:	Tag Free
Symbol:	GIDRP88
Synonyms:	C10orf28; GIDRP86; GIDRP88; PSORT
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for AF050198, the custom clone sequence may differ by one or more nucleotides

```

ATGCAGCAAGAATCAGAGAGATGCAGAGTTAGAGCCAGAAGGCCTGACATGGCACTTTATGTACCTAAAG
CTCGTAGGGGTGCAGTACTCCTTAAGACAGGTGATGAAGAAGAAAGCTGTGGTTCACCTAACTCTGTGGT
GAAAGAAAAGCAAAAAGAAAGTTCTCTCTCCCAAAAAGAAGTCTTTAAAGACAAACCGGAGGCTCGAAGA
CTAAATATCAATCCTGATAGAAAGGAGCATAATTGTAGAGAAGAAAAGAAATCTTCAACAAAATTAAGAA
TGGACACATGCCCTTCAAAAACAAATAGAGTTTGTCTAAGAGAGGAACCACTGAATCCAAAGAAGTATT
ATCCCAAGGACAACAGCAGGGAGCCCCAAATGCTGGGGTTATAACTAATGCACCTTTCAGAGACATTTT
AAACCAAGAAGGTGGAGTGTGGGAAGTTGAACTACGGATGTGACAGGACATGAGAGGATACTTCTTT
CACAGGCCTGTTTAGAAATCAGCGAGGCTCAAGTTCCAAGCAAACCAATCCAAAATGTGGAATTCTGTGA
CTTCAGTAGGCATGAACCTGATGGGGAAGCATTGGAAGACAAAGATTTGGAAGGCAGAATTGAAACTGAT
ACCAAGGTTTTGGAGATACTATATGAGTTTCCTAGAGTTTTAGTTCTGTGATGAAACCTGAGAATATGA
TTGTACCAATAAACTAAGCTCTGATTCTGAAATTGTACAACAAAGCATGCAAACATCAGATGGAATATT
GAATCCCAGCAGCGGAGGCATCACCCTACTTCTGTTCTGGAAGTCCAGATGGTGTCTTTGATCAAAC
TGCGTAGATTTTGAAGTTGAGAGTGTAGGTGGTATAGCCAATAGTACAGGTTTCATCTTAGATCAAAAGA
TACAGATTCCATTCCTGCAACTATGGGTCACATCTCTCTGTCAGAGAGCACAATGACACTGTTAGTCCA
GTAA
  
```


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for AF050198 unedited <pre> NGGGGTTCACAATTTGTATACGACTCACTATAGGCGGCCGCGAAATTCGCACGAGGGCTT CCCAAATGCTGGGATTACAGGTGAGGCTGCTGCAGAAAGACTTTTACACAGTTTGCCTTC AGAGACAGTCTATCGGCATGTTGTAGAGTTTTATTACTAAGAAAATAATGTTACTTACA TGATTGTGGTGGTCCATGCAGCAAGAATCAGAGAGATGCAGAGTTAGAGCCAGAAGGCC TGACATGGCACTTTATGTACCTAAAGCTCGTAGGGTGCAGTACTCCTTAAGACAGGTGA TGAAGAAGAAAGCTGTGGTTCACCTAACTCTGTGGTGAAAGAAAAGCAAAAAGAAAGTTC TCTCTCCCAAAAAGAAGTCTTTAAAGACAAACCGAGGCTCGAAGACTAAATATCAATCC TGATAGAAAGGAGCATAATTGTAGAGAAGAAAAGAAATCTTCAACAAAATTAAGAATGGA CACATGCCTTCAAAAACAAATAGAGTTTGTCTAAGAGAGGAACCACTGAATCCAAAGA AGTATTATCCCAAGGACAACAGCAGGGAGCCCCAAATGCTGGGGTTATACTAATGCACC TTTGCAGAGACATTTTAAACCAAGAAGGTGGAGTGTGGGAAGTTGAACTACGGATGT GACAGGACATGAGAGGATACTTCTTTCACAGGCCTGTTTAGAAATCAGCGAGGCTCAAGT TCCAAGCAAACCATTCAAAATGTTGGAATTCTGTGACTTCAGTAGGCATGAACCTGATG GGGGAGCATTTGAAGACAAAGATTGGAAGGCAGAATTGAACTGATACCCACGGTTTGG AGATACTATATGAGTTTCTAGAGTTTTAGTCCTGTGCATGAAACCTGAGATATGATTGT ACCATTAAGCTT </pre>
Restriction Sites:	NotI-NotI
ACCN:	AF050198
Insert Size:	3500 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	AF050198.1 , AAC05748.1
RefSeq Size:	1374 bp
RefSeq ORF:	1374 bp
Locus ID:	27291
Cytogenetics:	10q24.2