

Product datasheet for **SC105990**

ANKRD57 (SOWAHC) (AF085927) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKRD57 (SOWAHC) (AF085927) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKRD57
Synonyms:	HCCS1; hVps53L; PCH2E; pp13624
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for AF085927, the custom clone sequence may differ by one or more nucleotides

```
GAACGCCAGGAAGTCCATACTCTCAAACCCCTTTCTGCCCACTCACCTTCAGGTCCACAGCACTTTCCA
GCGTCATCTGAAATAGCAAAGGTGATGTGCAGCATCCCAGTGATGTGTTCTGTGTTATGAGAAGGAGTC
TGGGGTTCTGGATGCCTCAGGAACAGAATATGGACCTTAAGACCAGGCTGAGAGCCTGGCCTGGCACGC
AGCCTGCCCGCCCCGCGCCCTTCCTGTCCACCGCGGGCTCTCTGCAGCAGGGCTTGCCGAAGAGACTGC
TCCCGACCTCTCATGTCATCCTGCATCCAGAGAACGGCAGCATATCAAGCCGTTAACCTGATTCACTCG
GTTCCCTTGAAACCAGGCCTCAGTTAACCATGACCACAGCTACAGGCAAACCGCTGCTGTTGCTGCCTGC
TTTTGCTGAAGAATGGTATCTTGCCAGGCACGGCGGCTCATGCCTGTAGTCGCAACACTTTGGGAGACCA
AGGCAGGAGGATCACTTGAGTCCAGGAGTTTGAGACCAGCCTGGGCAACATACAGATGTCTACAAAATA
AAAATTTGCAAGGTGAAAAAAAAA
```



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for AF085927 unedited GCAGGAATAGGTAATACGATTTACTATAGGNNCGGCCGCGAAATTCGGCACGAGGAAAT ACAATTAATAAAAAAAAAAAGCAGGCAGGTGGNNCGGGTGCCTGTAATCTCAGCCTCC TGAGTATCGGATTACAGGCGTGAGCCACCGAGCCAGCCCTTCTCCATTTCTTTACTC CCACCTGCTCCTGGTTTCTAGGGAAGGGAACAGGAAGTTCACTTGGAAGACGTTTCC TAGGTTTTTCTCATCAGATGAAGCCATGGGAACGCCAGGAAGTCCATACTCTCAAACCCC TTTCTGCCCACTCACCTTCAGGTCCACAGCACTTTCAGCGTCATCTGAAATAGCAAAGG TGATGTGCAGCATCCCAGTGATGTGTTCTGTGTTATGAGAAGGAGTCTGGGGTTCTGGA TGCTCAGGAACAGAATATGGACCTTAAAGACCAGGCTGAGAGCCTGGCCTGGCACGCGAG CCTGCCGCCCGCCCTTCTGTCCACCGCGGGCTCTCTGCAGCAGGGCTTGCCGA AGAGACTGCTCCCGACCTCTCATGTCTCTGCATCCAGAGAACGGCAGCATATCAAGCC GTAAACCCTGATTCACTCGGTTCCCTTGAACACAGGCCTCAGTTAACCATGACCACAGCT ACAGGCAAACCGCTGCTGTTGCTGCCTGCTTTTGTGAAGATGATATCTTGCCANGCAC GGCGGGTCATGCCTGTAGTCGCAACACTTTGGGAGACCAAAGCACGANGGTACCCTGAG TTCACGAGTTTGAGACCAACCTGGGCAACATACAGATGTCTACAAAATTTAAATTTGCAA GGGGTGGTGGCATGCACCTGTAGTCCCACTGCTAGGGAAGCTGAGGTGGGCANAATGCTT GAACTCAGGAGTTTGAGACACGCTGGCCACCTGGCAAAACCACAATACAAAATGGCTG GTGGGGTGGCTCAG
Restriction Sites:	NotI-NotI
ACCN:	AF085927
Insert Size:	3000 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	AF085927.1
RefSeq Size:	585 bp
RefSeq ORF:	585 bp
Locus ID:	65124
Cytogenetics:	2q13