

Product datasheet for **SC106274**

FAM89A (BC009447) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM89A (BC009447) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM89A
Synonyms:	C1orf153; family with sequence similarity 89, member A; MGC15887; OTTHUMP00000036071; RP11-423F24.2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for BC009447, the custom clone sequence may differ by one or more nucleotides

ATGAGTGGGGCCCGGGCGGCGCCCGGGGCGGGCAACGGCGCGGTCCGGGGGCTGCGGGTGGACGGGC
 TGCCCCGCTGCCAAGAGCTTGAGCGGGCTGCTGCACTCGGCGTCGGGCGGCGGCGCTCTGGGGGCTG
 GCGGCACCTGGAGCGGCTGTACGCGCAGAAGTCGCGCATCCAGGACGAGCTGAGCCGCGGGGCCCCGGG
 GGCGGCGGGGCCCCGGGCGGCGAGCGCTGCCCGCAAGCCTCCCAACCTGGACGCCGCTCTGGCGCTGCTCC
 GCAAAGAGATGGTTGGTCTCCGCCAGCTGGACATGTCCTTGCTCTGCCAACTGTACAGCCTCTACGAGTC
 GATTTCAGGAGTACAAGGGGGCATGCCAGGCAGCCTCCAGCCAGACTGCACTTACGCTCTGGAGAACGGC
 TTCTTCGATGAAGAGGAGGAATATTTCCAGGAGCAGAAGTCCCTGCACGACAGGAGGGACCGAGGCCCTC
 CTCGGGACTTGCTACTGCCTGTCTCCTCCTCTCCAGCAGCGACTGGATTCTGGAGTCCATCTAG


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC009447 unedited <pre> NTGTTAGATTTGTAAACGACTCACTATAGGCGGCCGGAATCGGCACCAGGCGCGATGAG TGGGGCCCGGGCGGCGCCCGGGCCGCGGCAACGGCGCGGTCCGGGGGCTGCGGGTGGA CGGGCTGCCCCCGCTGCCAAGAGCTTGAGCGGGCTGCTGCACTCGGCGTCGGGCGGCGG CGCGTCTGGGGGCTGGCGGCACCTGGAGCGGCTGTACGCGCAGAAGTCGCGCATCCAGGA CGAGCTGAGCCGCGGGGGCCCGGGCGGCGGGGGCCGGGCGGCAGCGTGCCCGCCAA GCCTCCCAACCTGGACGCCGCTCTGGCGCTGCTCCGCAAAGAGATGGTTGGTCTCCGCCA GCTGGACATGTCTTGCTCTGCCAACTGTACAGCCTCTACGAGTCGATTACGGAGTACAA GGGGGCATGCCAGGCAGCTCCAGCCCAGACTGCACTTACGCTCTGGAGAACGGTCTCTT CGATGAAGAGGAGGAATATTTCCAGGAGCAGAAGTCCCTGCACGACAGGAGGGACCGAGG CCCTCCTCGGGACTTGCTACTGCCTGTCTCTCCCTCTCCAGCAGCGACTGGATTCTGGA GTCCATCTAGAGGGTCTTGGGAGGGATGTGACTGTTGGGAAGCCCTTCTACTGGACACG CTGTCATCATTTGCTGCTTCTTGCAAGAAAGCACCTCCGTTGTGGACGGTCTCGGGC ACAGGGGATGAGCGCTACCAAGTTTCATTTGTAGGCAGGGAGTTCTCCGCGGATGCATGGT GGCAGTCTGCTTTGATGGCAGCAGTTTCTGCTTAAGTGACCTAGAGGTCTCAGCAGTAT CCTCCACACCTATTTNATGNNAGTGCACTGCTGGGGATTTTCATATGAGAATATNACAGNA GGATCTCCGTGAAAGGCCTTANTGGGTGTTGTGTG </pre>
Restriction Sites:	NotI-NotI
ACCN:	BC009447
Insert Size:	1650 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC009447.1, AAH09447.1</u>
RefSeq Size:	1276 bp
RefSeq ORF:	555 bp
Locus ID:	375061
Cytogenetics:	1q42.2