

## Product datasheet for **SC107217**

### CAPON (NOS1AP) (NM\_014697) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CAPON (NOS1AP) (NM_014697) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAPON
Synonyms:	6330408P19Rik; CAPON; NPHS22
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_014697, the custom clone sequence may differ by one or more nucleotides

```

ATGCCTAGCAAAACCAAGTACAACCTTGTGGACGATGGGCACGACCTGCGGATCCCCTTGCAACAACGAGG
ACGCCTTCCAGCACGGCATCTGCTTTGAGGCCAAGTACGTAGGAAGCCTGGACGTGCCAAGGCCCAACAG
CAGGGTGGAGATCGTGGCTGCCATGCGCCGGATACGGTATGAGTTTAAAGCCAAGAACATCAAGAAGAAG
AAAGTGAGCATTATGGTTTCAGTGGATGGAGTGAAGTGATTCTGAAGAAGAAGAAAAGCTTCTTTTAT
TGCAGAAAAGGAATGGACGTGGGATGAGAGCAAGATGCTGGTATGCAGGACCCCATCTACAGGATCTT
CTATGTCTCTCATGATCCCAAGACTTGAAGATCTTCAGCTATATCGCTCGAGATGGTCCAGCAATATC
TTCAGGTGTAACGTCTTTAAATCCAAGAAGAAGAGCCAAGCTATGAGAATCGTTCGGACGGTGGGGCAGG
CCTTTGAGGTCTGCCACAAGCTGAGCCTGCAGCACACGCAGCAGAATGCAGATGGCCAGGAAGATGGAGA
GAGCGAGAGGAACAGCAACAGCTCAGGAGACCCAGGCCGCCAGCTCACTGGAGCCGAGAGGGCCTCCACG
GCCACTGCAGAGGAGACTGACATCGATGCGGTGGAGGTCCCACTTCCAGGGAATGATGTCCTGGAATTCA
GCCGAGGTGTGACTGATCTAGATGCTGTAGGGAAGGAAGGAGGCTCTCACACAGGCTCCAAGTTTCGCA
CCCCAGGAGCCCATGCTGACAGCCTCACCCAGGATGCTGCTCCCTTCTTCTTCCGAAGCCTCCAGGC
CTGGGCACAGAGACCCGCTGTCCACTCACCCAGATGCAGCTCCTCCAGCAGCTCCTCCAGCAGCAGC
AGCAGCAGACACAAGTGCTGTGGCCAGGTACTTGTGAAGGACCAGTTGGCTGCTGAGGCTGCGGC
GCGGCTGGAGGCCAGGCTCGCGTGCATCAGCTTTTGTGCAGAACAAGGACATGCTCCAGCACATCTCC
CTGCTGGTCAAGCAGGTGCAAGAGCTGGAAGTGAAGCTGTGAGGACAGAACGCCATGGGCTCCAGGACA
GCTTGCTGGAGATCACCTTCCGCTCCGGAGCCCTGCCCGTGTCTGTGACCCACGACCCCTAAGCCAGA
GGACCTGCATTCCGCCCGCTGGGCGGGGCTTGCTGACTTTGCCACCCCTGCGGGCAGCCCTTAGGT
AGGCGGACTGCTTGGTGAAGCTGGAGTGTCTTTCGCTTCTTCCGCCGAGGACACCCCGCCCCAGCGC
AGGGCGAGGCGCTCCTGGGCGGTCTGGAGCTCATCAAGTCCGAGAGTCAGGCATCGCTCGGAGTACGA
GTCCAACCGGACGAGAGCGAGGAGCGGACTCGTGGTCCAGGAGGAGCTGCCGCGCTGCTGAATGTC
CTGCAGAGGCAGGAAGTGGGCGACGGCCTGGATGATGAGATCGCCGTGTAG

```



[View online >](#)

<b>Restriction Sites:</b>	Please inquire
<b>ACCN:</b>	NM_014697
<b>Insert Size:</b>	2600 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_014697.1</a> , <a href="#">NP_055512.1</a>
<b>RefSeq Size:</b>	2891 bp
<b>RefSeq ORF:</b>	1521 bp
<b>Locus ID:</b>	9722
<b>UniProt ID:</b>	<a href="#">O75052</a>
<b>Cytogenetics:</b>	1q23.3
<b>Gene Summary:</b>	<p>This gene encodes a cytosolic protein that binds to the signaling molecule, neuronal nitric oxide synthase (nNOS). This protein has a C-terminal PDZ-binding domain that mediates interactions with nNOS and an N-terminal phosphotyrosine binding (PTB) domain that binds to the small monomeric G protein, Dexas1. Studies of the related mouse and rat proteins have shown that this protein functions as an adapter protein linking nNOS to specific targets, such as Dexas1 and the synapsins. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2009]</p> <p>Transcript Variant: This variant (1) encodes the full-length protein (isoform 1). Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>