

5' Read Nucleotide Sequence:	>OriGene 5' read for NM_198404 unedited <pre> TTGTAATACGACTCACTATAGGGCGGCCGCAATTCGGCACGAGGGGCTGGGGCTGAATA AAGGCGTGCTTGTGTGGTAGTGTCTCTTTTAAAAATCTCAAAGCCAAGAACAAGC TGAAATAGCATCTTCAAAAAATGGAGCGTAAATAAACAGAGAGAAAAAGAAAGGAGT ATGAAGGGAAACACAACAGCCTGGAAGATACTGATCAAGGAAAGAACTGCAATCCACAC TGATGACCCTCAACGTTGGTGGATATTTATACATTACTCAAAAACAACTGACCAAGT ACCCAGACACTTTCCTTGAAGGTATAGTAAATGGAAAAATCCTCTGCCCGTTTGATGCTG ATGGTCATTATTTATAGACAGGGATGGTCTCCTCTTCAGGCATGTCTAAACTTCCTAC GAAATGGAGAAGTCTATTGCCCGAAGGGTTTCGAGAAATCAACTTCTTGCACAAGAAG CAGAAATCTTTCAGCTCAAGGGACTGGCAGAGGAAGTGAAATCCAGTGGGAGAAAGAAC AGCTAACACCCAGAGAGACTACTTCTTGGAAATAACAGATAACCACGATCGTTCACAAG GATTAAGAATCTTCTGTAATGCTCCTGATTTATATCAAAAATAAAGTCTCGCATTGTTT TGGTGTCCAAAAGCAGGCTGGATGGATTCCAGAGGAGTTTCAATATCGTCAAATATCA TCCAATTTAAATACTTCATAAAGTCTGAAAATGGCACTCGACTTGTACTAAAGGAAAGAC ACACCTTTGTCTGTACCTTGGNAACTCTAGTTTGAGGCTATCATGATGGCTTTAAAGT GTGGCTTTAACTGCTGACCAGCCTGNNATGTTCCAAAGGGTCAATNGTCCAG </pre>
Restriction Sites:	NotI-NotI
ACCN:	NM_198404
Insert Size:	4000 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198404.1 , NP_940686.1
RefSeq Size:	1878 bp
RefSeq ORF:	780 bp
Locus ID:	386618
UniProt ID:	Q8WVF5
Cytogenetics:	13q14.12-q14.13