

Product datasheet for **SC108409**

Cullin 1 (CUL1) (NM_003592) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cullin 1 (CUL1) (NM_003592) Human Untagged Clone
Tag:	Tag Free
Symbol:	Cullin 1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_003592, the custom clone sequence may differ by one or more nucleotides

```
ATGTCGTC A ACCCGAGCCAGA ACCCCACGGCCTGAAGCAGATTGGCCTGGACCAGATCTGGGACGACC
TCAGAGCCGGCATCCAGCAGGTGTACACACGGCAGAGCATGGCCAAGTCCAGATATATGGAGCTCTACAC
TCATGTTTATAACTACTGTACTAGTGTTCACCAGTCAAACCAAGCACGAGGAGCTGGAGTTCCTCCTTCT
AAGTCGAAAAAGGGGCAGACACCTGGAGGAGCTCAGTTTGTGGCCTGGAATTATATAAACGACTTAAGG
AATTTTTGAAGAATTACTTGACAATCTCTTAAGGATGGAGAAGATTTGATGGATGAGAGTGTACTGAA
ATTCTACACTCAACAATGGGAAGATTATCGATTTTCAAGCAAAGTGTGAATGGAATTTGTGCCTACCTC
AATAGACATTGGGTTCCGCGTGAATGTGACGAAGGACGAAAAGGAATATATGAAATCTATTGCTTGCAT
TGGTGACTTGGAGAGACTGTCTGTTGAGGCCACTGAATAAACAGGTAACAAATGCTGTTTTAAAGCTGAT
TGAAAAGGAAAGGAATGGTGAACCATCAATACAAGATTGATTAGTGGAGTTGTACAGTCTTACGTGGAA
TTGGGGCTGAATGAAGATGATGCATTTGCAAAGGGCCCTACGTTAACAGTGTATAAAGAATCCTTTGAAT
CTCAATTTTTGGCTGACACAGAGAGATTTATACCAGAGAGAGTACTGAATTCCTGCAGCAGAACCAGT
TACTGAATATATGAAAAAGGCAGAGGCTCGTCTGCTTGAGGAACAACGAAGATTCAGGTTTACCTTCAT
GAAAGCACACAAGATGAATTAGCAAGGAAATGTGAACAAGTCTCATTGAAAAACTTGGAAATTTTCC
ACACAGAAATTCAGAATTTATTGGATGCTGACAAAAATGAAGATTTGGGACGCATGTATAATCTTGATC
TAGAATCCAGGATGGCCTAGGAGAATTGAAAAACTGTTGGAGACACACATTCATAATCAGGGTCTTGCA
GCCATTGAAAAGTGTGGAGAAGCTGCTTAAATGACCCCAAAATGTATGTACAGACAGTCTTGATGTTT
ATAAAAAACAATGCCCTGGTAATGTCTGCATTCAACAATGACGCTGGCTTTGTGGCTGCTCTTGATAA
GGCTTGTGGTTCGCTTATAAACAACAACGCGGTTACCAAGATGGCCCAATCATCCAGTAAATCCCCTGAG
TTGCTGGCTCGATACTGTGACTCCTTGTGAAGAAAAGTTCCAAGAACCAGAGGAGGCAGAACTAGAAG
ACACACTCAATCAAGTGATGGTTGCTTCAAGTACATAGAAGACAAAGAGTATTTGAGAAGTCTATGC
GAAGATGCTCGCCAAGAGGCTCGTCCACAGAACAGTGAAGTACGATGCCGAAGCCAGCATGATCTCC
AAGTTAAAGCAAGCTTGCGGGTTGAGTACACCTCTAACTTCAGCGCATGTTTCAAGACATTGGCGTGA
GCAAAGATCTGAACGAGCAATTCAAAAAGCACTTGACAAACTCAGAACCCTAGACTTGGATTTGAGCAT
TCAAGTGTGAGTCCGGGTCCTGGCCCTTCCAGCAGTCTGTACATTTGCCTTGGCGTCAGAGTTGGAA
CGTAGTTATCAGCGATTCACAGCTTTCTACGCCAGCCGCCACAGTGGCCGAAAATGACGTGGTTATATC
AGTTGTCTAAAGGAGAATTGGTAACTAAGTCTTCAAAAAAGATATACTTTGCAGGCGTCGACATTCCA
GATGGCTATCCTGCTTCAGTACAACACGGAAGATGCCTACACTGTGCAGCAGCTGACCGACAGCACTCAA
ATTAATAATGGACATTTTGGCGCAAGTTTTACAGATTTTATTAAGTCAAGCTATTGGTCTTGAAGATG
AAAATGCAATGTTGATGAGGTGGAATTAAGCCAGATACCTTAATAAAATATATCTTGTTTATAAAAA
TAAGAAATTAAGGGTTAACATCAATGTGCCAATGAAAACCGAACAGAAGCAGGAACAAGAAACCACACAC
AAAAACATCGAGGAAGACCGCAAACTACTGATTGAGGCGCCATCGTGAGAATCATGAAGATGAGGAAGG
TTCTGAAACACCAGCAGTTACTTGGCGAGGTCCTCACTCAGCTGTCTCCAGGTTCAAACCTCGAGTCCC
TGTGATCAAGAAATGCATTGACATTCTAATTGAGAAAAGATTTGGAGCGAGTGGATGGTGAAGGAC
ACCTACAGTACTTGGCTTAA
```

5' Read Nucleotide Sequence:	<p>>OriGene 5' read for NM_003592 unedited NAATTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGGCACGAGGCGGAAGTCCGC GCCGATTAGGAGGGCCACCTCTGACTTTCTGATCTGGGTGCAACCCCACTCCCTGTGC CGGTGACAAAGGTTTGCCTGCAATGAGATTTCTTCTACATTTAAAGGACATCCTTT CTGAGCTGCTGTGAATAAATTTGGAATGGTACTGTATATTTTCATCTAATGGAGAAGTAC CTGTACTTTGAATAAGGATTGCTGCACTGGACGACTTTAGAACATCCCTCACAATGTCGT CAACCCGGAGCCAGAACCCCAAGGCTGAAGCAGATTGGCCTGGACCAGATCTGGGACG ACCTCAGAGCCGCATCCAGCAGGTGTACACACGGCAGAGCATGGCCAAGTCCAGATATA TGGAGCTCTACACTCATGTTTATAACTACTGTAAGTGTTCACCAAGTCAAACCAAGCAC GAGGAGCTGGNAGTTCCTCCTTCTAAGTCGAAAAAGGGGACACACCTGNAGGAGCTCAG TTTGTTGGCCTGGAATTATATAACGACTTAAGGAATTTTGAAGAATTACTGNACAAT CTTCTTAANGATGGAGAAGATTTGATGATGAAAAGTGTACCTGAAATTCTACTNACAC ATGGGAAGATANATCGATNTCAAACAAAGTGTGAATGGAAATTTGTGCCTACTCATAAC AAATTGGNTTTCGCCGTGATGTNNACGAAGGACGAAAGGATATNTGAAATCTATTCT TGCATTGGGACTGNAAAAAATCTGNTCTTGTTCAGCACTGAATACAGGCAAAGGCCTG</p>
3' Read Nucleotide Sequence:	<p>>OriGene 3' read for NM_003592 unedited CTTGAACGCGCCGAATCTAGGATCGAGTTTTTTTTTTTTTTTTTTTGCACAATCATTAT ATCTTTTATTTACATATAGAAATTAGTAAACTTAAAACAAAACCAATGCGACTGTAC AGAAATCCTTTAATTTTCTTTATTTTACACATTAATAAATGAAACACACTATACAAATG GTTTTCATAGCAGTTACATGTGGTCCACTCAGACTTGCTCTCAAAGTGTGCCAACGT GACAAGTGTGCCAAGACGACAGCTGCGAATCGGGGAGCCATGTATTGACTTAGTGCGG TAGCATGCAGCCCTCGGCACGGAGGGTCACTCCCGCTGTTCTTAAATTAGGGTAAATGG CGTTGGTGTCCGATTTACAGACTTGAAATGTATCAATCCTGAGCAGTCCGATGTAAT CCGAGACCAGCTGATGGGAGTCTTGGTCTCAGTTTTTAAAGGGAGGTCCAATGGCGGCA GGCTGGTGTATGAACTTGAGTTGTTTTTATTCTTTCCAACATGAACTATTTGTGCGG GTCACACAGTCAGACCCTCCAGAAGGGTTAAGCCAAGTAACTGTAGGTGCTTTTTAC CATCCACTCGTCCAAATATTCTTTCTCAATTAGAATGTCAATGCATTTCTTGATCAG GGACTCGAGGTTTGAACCTGGAGGACAGCTGAGTGAGGACCTCGCAAGTAACTGCTGGT GTTTCAGAACCTTCCTCATCTTCAATGATTCTCACGATGGCCGCTGAATCAGTAGTTG GGGCCTCCTCGATGTTNTTGTGTGTTTCTTGGTCTGCTTCTGNNTCGNTNTCATTG CACATTGATGNTAACCTTAATTNCTTATTTTATAACCAAGAATAATNNTATAAGGTATC TGGCTTATTCCACTCATCATTGCATTTCTTCAAGACCATAGCTCACTTTATAAA TCTGTAAGTGGCAAAGTCAATNNATTGGAGG</p>
Restriction Sites:	NotI-NotI
ACCN:	NM_003592
Insert Size:	2840 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003592.2 , NP_003583.2
RefSeq Size:	3226 bp
RefSeq ORF:	2331 bp
Locus ID:	8454
UniProt ID:	Q13616
Cytogenetics:	7q36.1
Domains:	CULLIN
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Oocyte meiosis, TGF-beta signaling pathway, Ubiquitin mediated proteolysis, Wnt signaling pathway

Gene Summary:

Core component of multiple cullin-RING-based SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complexes, which mediate the ubiquitination of proteins involved in cell cycle progression, signal transduction and transcription. SCF complexes and ARIH1 collaborate in tandem to mediate ubiquitination of target proteins (PubMed:27565346). In the SCF complex, serves as a rigid scaffold that organizes the SKP1-F-box protein and RBX1 subunits. May contribute to catalysis through positioning of the substrate and the ubiquitin-conjugating enzyme. The E3 ubiquitin-protein ligase activity of the complex is dependent on the neddylation of the cullin subunit and exchange of the substrate recognition component is mediated by TIP120A/CAND1. The functional specificity of the SCF complex depends on the F-box protein as substrate recognition component. SCF(BTRC) and SCF(FBXW11) direct ubiquitination of CTNNB1 and participate in Wnt signaling. SCF(FBXW11) directs ubiquitination of phosphorylated NFKBIA. SCF(BTRC) directs ubiquitination of NFKBIB, NFKBIE, ATF4, SMAD3, SMAD4, CDC25A, FBXO5 and probably NFKB2. SCF(BTRC) and/or SCF(FBXW11) direct ubiquitination of CEP68 (PubMed:25704143, PubMed:25503564). SCF(SKP2) directs ubiquitination of phosphorylated CDKN1B/p27kip and is involved in regulation of G1/S transition. SCF(SKP2) directs ubiquitination of ORC1, CDT1, RBL2, ELF4, CDKN1A, RAG2, FOXO1A, and probably MYC and TAL1. SCF(FBXW7) directs ubiquitination of CCNE1, NOTCH1 released notch intracellular domain (NICD), and probably PSEN1. SCF(FBXW2) directs ubiquitination of GCM1. SCF(FBXO32) directs ubiquitination of MYOD1. SCF(FBXO7) directs ubiquitination of BIRC2 and DLGAP5. SCF(FBXO33) directs ubiquitination of YBX1. SCF(FBXO1) directs ubiquitination of BCL6 and DTL but does not seem to direct ubiquitination of TP53. SCF(BTRC) mediates the ubiquitination of NFKBIA at 'Lys-21' and 'Lys-22'; the degradation frees the associated NFKB1-RELA dimer to translocate into the nucleus and to activate transcription. SCF(CCNF) directs ubiquitination of CCP110. SCF(FBXL3) and SCF(FBXL21) direct ubiquitination of CRY1 and CRY2. SCF(FBXO9) directs ubiquitination of TT11 and TELO2. SCF(FBXO10) directs ubiquitination of BCL2.[UniProtKB/Swiss-Prot Function]