

Product datasheet for **SC108712**

FKBP14 (NM_017946) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FKBP14 (NM_017946) Human Untagged Clone
Tag:	Tag Free
Symbol:	FKBP14
Synonyms:	EDSKMH; EDSKSCL2; FKBP22; IPBP12
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF within SC108712 sequence for NM_017946 edited (data generated by NextGen Sequencing) ATGAGGCTTTTCTGTGGAACGCGGTCTTGACTCTGTTCTGCACTTCTTTGATTGGGGCT TTGATCCCTGAACCAGAAGTGAAAATTGAAGTTCTCCAGAAGCCATTCATCTGCCATCGC AAGACCAAAGGAGGGGATTTGATGTTGGTCCACTATGAAGGCTACTTAGAAAAGGACGGC TCCTTATTTCACTCCACTCACAACATAACAATGGTCAGCCATTTGGTTTACCCTGGGC ATCCTGGAGGCTCTCAAAGGTTGGGACCAGGGCTTGAAGGAATGTGTGTAGGAGAGAAG AGAAAGCTCATCATTCTCCTGCTCTGGGCTATGGAAAAGAAGGAAAAGTAAAATTCCC CCAGAAAGTACTGATATTTAATATTGATCTCCTGGAGATTCGAAATGGACCAAGATCC CATGAATCATTCCAAGAAATGGATCTTAATGATGACTGGAACTCTCTAAAGATGAGGTT AAAGCATATTTAAGAAGGAGTTTGAAAAACATGGTGCGGTGGTGAATGAAAGTCATCAT GATGCTTTGGTGGAGGATATTTTGATAAAGAAGATGAAGACAAAGATGGGTTTATATCT GCCAGAGAATTTACATATAAACACGATGAGTTATAG Clone variation with respect to NM_017946.2



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017946.2](#), [NP_060416.1](#)

RefSeq Size: 2247 bp

RefSeq ORF: 636 bp

Locus ID: 55033

UniProt ID: [Q9NWM8](#)

Cytogenetics: 7p14.3

Domains: FKBP, EFh

Protein Families: Druggable Genome

Gene Summary: The protein encoded by this gene is a member of the FK506-binding protein family of peptidyl-prolyl cis-trans isomerases. The encoded protein is found in the lumen of the endoplasmic reticulum, where it is thought to accelerate protein folding. Defects in this gene are a cause of a type of Ehlers-Danlos syndrome (EDS). Both a protein-coding variant and noncoding variants are transcribed from this gene. [provided by RefSeq, Mar 2012]
Transcript Variant: This variant (1) represents the protein-coding transcript. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.