

Product datasheet for SC108827

HIP1 (NM_005338) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HIP1 (NM_005338) Human Untagged Clone
Tag:	Tag Free
Symbol:	HIP1
Synonyms:	HIP-I; ILWEQ; SHON; SHONbeta; SHONGgamma
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF within SC108827 sequence for NM_005338 edited (data generated by NextGen Sequencing)

```

ATGGATCGGATGGCCAGCTCCATGAAGCAGGTGCCCAACCCACTGCCCAAGGTGCTGAGC
CGGCGCGGGTTCGGCGCTGGGCTGGAGGCGGCGAGCGGAGAGCTTCGAGCGGACTCAG
ACTGTGTCAGCATCAATAAGGCCATTAATACGCAGGAAGTGCTGTAAGGAAAAACACGCC
AGAACGTGCATACTGGGCACCCACCATGAGAAAGGGGCACAGACCTTCTGGTCTGTTGTC
AACCGCTGCCTCTGTCTAGCAACGCAGTGCTCTGCTGGAAGTTCTGCCATGTGTTCCAC
AAACTCCTCCGAGATGGACACCCGAACGTCCTGAAGGACTCTCTGAGATACAGAAATGAA
TTGAGTGACATGAGCAGGATGTGGGGCCACCTGAGCGAGGGGTATGGCCAGCTGTGCAGC
ATCTACCTGAAACTGCTAAGAACCAAGATGGAGTACCACACCAAAAATCCCAGGTTCCCA
GGCAACCTGCAGATGAGTGACCGCCAGCTGGACGAGGCTGGAGAAAGTGACGTGAACAAC
TTTTTCCAGTTAACAGTGGAGATGTTTGACTACCTGGAGTGTGAACTCAACCTCTTCCAA
ACAGTATTCAACTCCCTGGACATGTCCCCTCTGTGTCCGTGACGGCAGCAGGGCAGTGC
CGCCTCGCCCCGCTGATCCAGGTCATCTTGGACTGCAGCCACCTTTATGACTACACTGTC
AAGCTTCTCTTCAAACCTCCACTCCTGCCTCCCAGCTGACACCTGCAAGGCCACCGGGAC
CGCTTCATGGAGCAGTTTACAAAGTTGAAAGATCTGTTCTACCGTCCAGCAACCTGCAG
TACTTCAAGCGGCTCATTAGATCCCCAGCTGCCTGAGAACCCACCAACTTCTGCGA
GCCTCAGCCCTGTCAGAACATATCAGCCCTGTGGTGGTGATCCCTGCAGAGGCCTCATCC
CCCGACAGCGAGCCAGTCTAGAGAAGGATGACCTCATGGACATGGATGCCTCTCAGCAG
AATTTATTTGACAACAAGTTTGATGACATCTTTGGCAGTTTATTAGCAGTGTATCCCTTC
AATTTCAACAGTCAAATGGTGTGAACAAGGATGAGAAGGACCCTTAATTGAGCGACTA
TACAGAGAGATCAGTGGATTGAAGGCACAGCTAGAAAACATGAAGACTGAGAGCCAGCGG
GTTGTGCTGCAGCTGAAGGGCCACGTCAGCGAGCTGGAAGCAGATCTGGCCGAGCAGCAG
CACCTGCGGCAGCAGGCGGCCGACGACTGTGAATTCTGCGGGCAGAACTGGACGAGCTC
AGGAGGCAGCGGGAGGACACCGAGAAGGCTCAGCGGAGCCTGTCTGAGATAGAAAGGAAA
GCTCAAGCCAATGAACAGCGATATAGCAAGCTAAAGGAGAAGTACAGCGAGCTGGTTCAG
AACCACGCTGACCTGCTGCGGAAGAATGCAGAGGTGACCAACAGGTGTCCATGGCCAGA

```



[View online »](#)

```

CAAGCCCAGGTAGATTTGGAACGAGAGAAAAAGAGCTGGAGGATTTCGTTGGAGCGCATC
AGTGACCAGGGCCAGCGGAAGACTCAAGAACAGCTGGAAGTTCTAGAGAGCTTGAAGCAG
GAACTTGGCCACAAGCCAACGGGAGCTTCAGGTTCTGCAAGGCAGCCTGGAAACTTCTGCC
CAGTCAGAAGCAAACCTGGGCAGCCGAGTTCGCCGAGCTAGAGAAGGAGCGGGACAGCCTG
GTGAGTGGCGCAGCTCATAGGGAGGAGGAATTATCTGCTCTTCGGAAAGAACTGCAGGAC
ACTCAGCTCAAACCTGGCCAGCACAGAGGAATCTATGTGCCAGCTTGCCAAAGACCAACGA
AAAATGCTTCTGGTGGGTCCAGGAAGGCTGCGGAGCAGGTGATACAAGACGCCCTGAAC
CAGCTTGAAGAACCCTCCTCTCATCAGCTGCGCTGGGTCTGCAGATCACCTCCTCTCCACG
GTCACATCCATTTCCAGCTGCATCGAGCAACTGGAGAAAAGCTGGAGCCAGTATCTGGCC
TGCCCAAGACATCAGTGGACTTCTCCATTCCATAAACCCTGCTGGCCACTTGACCAGC
GACGCCATTGCTCATGGTGCCACCCTGCCTCAGAGCCCCACCTGAGCCTGCCGACTCA
CTGACCGAGGCCTGTAAGCAGTATGGCAGGAAACCTCGCTACCTGGCCTCCCTGGAG
GAAGAGGGAAGCCTTGAGAATGCCGACAGCACAGCCATGAGGAACTGCCTGAGCAAGATC
AAGGCCATCGGCGAGGAGCTCCTGCCAGGGGACTGGACATCAAGCAGGAGGAGCTGGGG
GACCTGGTGGACAAGGAGATGGCGGCCACTTCAGCTGCTATTGAAACTGCCACGGCCAGA
ATAGAGGAGATGCTCAGCAAATCCCAGAGAGGAGACACAGGAGTCAAATTGGAGGTGAAT
GAAAGGATCCTTGGTTGCTGTACCAGCCTCATGCAAGCTATTCAGGTGCTCATCGTGCC
TCTAAGGACCTCCAGAGAGAGATTGTGGAGAGCGGCAGGGGTACAGCATCCCCTAAAGAG
TTTTATGCCAAGAACTCTCGATGGACAGAAGGACTTATCTCAGCCTCCAAGGCTGTGGGC
TGGGGAGCCACTGTATGGTGGATGCAGCTGATCTGGTGGTACAAGGCAGAGGGAATTT
GAGGAGCTAATGGTGTGTCTCATGAAATTGCTGCTAGCACAGCCAGCTTGTGGCTGCA
TCCAAGGTGAAAGCTGATAAGGACAGCCCAACCTAGCCAGCTGCAGCAGGCCTCTCGG
GGAGTGAACCAAGCCACTGCCGGCTTGTGGCCTCAACCATTTCCGGCAAATCACAGATC
GAAGAGACAGACAACATGGACTTCTCAAGCATGACGCTGACACAGATCAAACGCCAAGAG
ATGGATTCTCAGTTAGGGTCTAGAGCTAGAAAATGAATTGCAGAAGGAGCGTCAAAAA
CTGGGAGAGCTTCGAAAAAGCACTACGAGCTTCTGGTGGTCTGAGGGCTGGGAAGAA
GGAACAGAGGCATCTCCACCTACTGCAAGAAGTGGTAACCGAAAAAGAATAG

```

Clone variation with respect to NM_005338.5

Restriction Sites:

NotI-NotI

ACCN:

NM_005338

Insert Size:

7000 bp

OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_005338.3](#), [NP_005329.2](#)

RefSeq Size: 6626 bp

RefSeq ORF: 2745 bp

Locus ID: 3092

UniProt ID: [O00291](#)

Cytogenetics: 7q11.23

Domains: I_LWEQ

Protein Families: Druggable Genome

Protein Pathways: Huntington's disease

Gene Summary: The product of this gene is a membrane-associated protein that functions in clathrin-mediated endocytosis and protein trafficking within the cell. The encoded protein binds to the huntingtin protein in the brain; this interaction is lost in Huntington's disease. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.