

# Product datasheet for SC108966

## ATP6V1G2 (NM\_138282) Human Untagged Clone

### **Product data:**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type	Expression Plasmids
Dreduct Norses	ATRC//1C2 (NM 128282) Using Untergred Class
Product Name:	ATPOVIGE (NM_138282) Human Unlagged Clone
Tag:	Tag Free
Symbol:	ATP6V1G2
Synonyms:	ATP6G; ATP6G2; NG38; VMA10
Mammalian Cell Selection:	None
Vector:	pCMV6-XL5
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	<pre>&gt;OriGene ORF within SC108966 sequence for NM_138282 edited (data generated by NextGen Sequencing) ATGGAGGTGGAGCAATACCGCAGAGAGCGAGAGCACGAATTCCAGAGCAAGCA</pre>
5' Read Nucleotide Sequence:	<pre>&gt;DriGene 5' read for NM_138282 unedited GTAATACGACTCACTTATAGGGCGGCCGCGAATTCGCACGAGGGGCAGCTGAGAAGGTGG CAGATGCCAGAAAGAGGAAGGCCCGGCGACTGAAGCAGGCAG</pre>
Restriction Sites:	Notl-Notl



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	TP6V1G2 (NM_138282) Human Untagged Clone – SC108966
ACCN:	NM_138282
Insert Size:	4000 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Met	<ul> <li>hod: 1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ul>
RefSeq:	<u>NM 138282.1, NP 612139.1</u>
RefSeq Size:	1440 bp
RefSeq ORF:	234 bp
Locus ID:	534
UniProt ID:	<u>095670</u>
Cytogenetics:	6p21.33
Domains:	V-ATPase_G
Protein Pathways:	Epithelial cell signaling in Helicobacter pylori infection, Metabolic pathways, Oxidative phosphorylation, Vibrio cholerae infection

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

#### SC108966 ATP6V1G2 (NM\_138282) Human Untagged Clone – SC108966

Gene Summary: This gene encodes a component of vacuolar ATPase (V-ATPase), a multisubunit enzyme that mediates acidification of intracellular compartments of eukaryotic cells. V-ATPase dependent acidification is necessary for such intracellular processes as protein sorting, zymogen activation, receptor-mediated endocytosis, and synaptic vesicle proton gradient generation. V-ATPase is composed of a cytosolic V1 domain and a transmembrane V0 domain. The V1 domain consists of three A and three B subunits, two G subunits plus the C, D, E, F, and H subunits. The V1 domain contains the ATP catalytic site. The V0 domain consists of five different subunits: a, c, c', c'', and d. Additional isoforms of many of the V1 and V0 subunit proteins are encoded by multiple genes or alternatively spliced transcript variants. This encoded protein is one of three V1 domain G subunit proteins. This gene had previous gene symbols of ATP6G and ATP6G2. Alternatively spliced transcript variants encoding different isoforms have been described. Read-through transcription also exists between this gene and the downstream DEAD (Asp-Glu-Ala-Asp) box polypeptide 39B (DDX39B) gene. [provided by RefSeq, Feb 2011]

Transcript Variant: This variant (2) lacks a segment of the 5' UTR and 5' coding region, and uses a downstream start codon, compared to variant 1. The encoded isoform (b), also known as ATP6Galt, is shorter at the N-terminus, compared to isoform a.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US