

Product datasheet for **SC109269**

C6orf62 (NM_030939) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C6orf62 (NM_030939) Human Untagged Clone
Tag:	Tag Free
Symbol:	C6orf62
Synonyms:	dj30M3.2; Nbla00237; XTP12
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF within SC109269 sequence for NM_030939 edited (data generated by NextGen Sequencing)

```
ATGGGGACCCAACTCCCGAAGAAACAAGCTCTGAACAGACTACGTGCTCAGCTTAGA
AAGAAAAAGAATCTCTAGCTGACCAGTTTGACTTCAAGATGTATATTGCCTTTGTATTC
AAGGAGAAGAAGAAAAGTCAGCACTTTTTGAAGTGTCTGAGGTTATACCAAGTCATGACA
AATAATTATGAAGAAATATCCTGAAAGGTGTGCGAGATTCCAGCTATTCCTTGGAAAGT
TCCCTAGAGCTTTTACAGAAGGATGTGGTACAGCTCCATGCTCCTCGATATCAGTCTATG
AGAAGGGATGTAATTGGCTGTAAGGAGATGGATTTTCATTCTTTGGCCTCGGAATGAT
ATTGAAAAATCGTCTGTCTCCTGTTTTCTAGGTGAAAGAATCTGATGAGCCTTTTAGG
CCTGTTACAGCCAAATTTGAGTTTCATCATGGTACTATGAAAAACAGTTTCTGCATGTA
CTGAGCCGCAAGGACAAGACTGGAATCGTTGTCAACAATCCTAACCAGTCAGTGTCTCTC
TTCATTGACAGACAGCACTTGCAGACTCCAAAAACAAAAGCTACAATCTTCAAGTTATGC
AGCATCTGCCTCTACCTGCCACAGGAACAGCTCACCCACTGGGCAGTTGGCACCATAGAG
GATCACCTCCGTCCTTATATGCCAGAGTAG
```

Clone variation with respect to NM_030939.4



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030939.3](#), [NP_112201.1](#)

RefSeq Size: 2424 bp

RefSeq ORF: 690 bp

Locus ID: 81688

UniProt ID: [Q9GZU0](#)

Cytogenetics: 6p22.3