

Product datasheet for **SC110009**

CPNE1 (NM_152927) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CPNE1 (NM_152927) Human Untagged Clone
Tag:	Tag Free
Symbol:	CPNE1
Synonyms:	COPN1; CPN1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_152927, the custom clone sequence may differ by one or more nucleotides

```

ATGGCCACTGCGTGACCTTGGTTACAGCTGTCCATTTCTGTGACCATCTCATTGACAAGGACATCGGCT
CCAAGTCTGACCCACTCTGCGTCCTTTTACAGGATGTGGGAGGGGGCAGCTGGGCTGAGCTTGGCCGGAC
TGAACGGGTGCGGAAGTCTCAAGCCCTGAGTTCTCAAGACTCTACAGCTTGAGTACCGCTTTGAGACA
GTCCAGAAGCTACGCTTTGGAATCTATGACATAGACAACAAGACGCCAGAGCTGAGGGATGATGACTTCC
TAGGGGGTGTGAGTGTTCCCTAGGACAGATTGTGCCAGCCAGGTAAGTACTGACTCTCCCCTTGATGCTGAA
GCCTGGAAAACCTGCTGGGCGGGGGACCATCACGGTCTCAGCTCAGGAATTAAGGACAATCGTGTAGTA
ACCATGGAGGTAGAGGCCAGAACTAGATAAGAAGGACTTCTGGGAAAATCAGATCCATTTCTGGAGT
TCTTCCGCCAGGGTGTGGGAAATGGCACCTGGTGTACAGATCTGAGGTCATCAAGAACAACCTGAACCC
TACATGGAAGCGTTTCTCAGTCCCCGTTACGATTTCTGTGGTGGGAACCCAGCACACCCATCCAGGTG
CAATGCTCCGATTATGACAGTGACGGGTACATGATCTCATCGGTACCTCCACACCAGCTTGGCCAGC
TGCAGGCAGTCCCGGCTGAGTTTGAATGCATCCACCCTGAGAAGCAGCAGAAAAAGAAAGCTACAAGAA
CTCTGGAAGTATCCGTGTCAAGATTTGTGGGTAGAAAACAGAGTACTCCTTTCTGGACTATGTGATGGGA
GGCTGTACAGTCAACTTCACTGTGGGCGTGGACTTCACTGGCTCCAATGGAGACCCCTCCTCACCTGACT
CCCTACACTACCTGAGTCCAACAGGGGTCAATGAGTACCTGATGGCACTGTGGAGTGTGGGCAGCGTGGT
TCAGGACTATGACTCAGACAAGCTGTCCCTGCATTTGGATTTGGGGCCAGGTTCCCTCTGACTGGCAG
GTCTCGCATGAATTTGCCTTGAATTTCAACCCAGTAACCCCTACTGTGCAGGCATCCAGGGCATTGTGG
ATGCCTACCGCAAGCCCTGCCCAAGTTCGCCTCTATGGCCCTACCAACTTTGCACCCATCATCAACCA
TGTGGCCAGGTTGCAGCCAGGCTGCACATCAGGGGACTGCCTCGCAATACTTCATGCTGTTGCTGCTG
ACTGATGGTGTGTGACGGATGTGGAAGCCACACGTGAGGCTGTGGTGGTGCCTCGAAGCTGCCATG
CAGTGATCATTGTGGGTGTGGGTGGTGTGACTTTGAGGCCATGGAGCAGCTGGAGCTGATGGTGGACC
CCTGCATACACGTTCTGGCAGGCTGTGCCCGGACATTGTGCAGTTTGTACCTACCGCCGGTCCAG
AATGCCCTCGGGAGGCATTGGCACAGACCGTGTCTGCAGAAGTGCCACACAAGTGGTCTCATACTTCA
GGGCCAGGGTTGGGCCCGCTCAAGCCACTTCCACCCTCAGCCAAGGATCTGCACAGGCCCCCGAGCC
CTAG
    
```

5' Read Nucleotide Sequence:

>OriGene 5' read for NM_152927 unedited

```

GATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGGCACGAGGGCTGCGTTGGG
CTCCGGGAAGCCGTTTCGGGCTGGGGCTGTCCGGCCGCGGGGCGGAGGCACTCGCGCGGGGG
ATGGCCACTGCGTGACCTTGGTTACAGCTGTCCATTTCTGTGACCATCTCATTGACAAG
GACATCGGCTCCAAGTCTGACCCACTCTGCGTCCTTTTACAGGATGTGGGAGGGGGCAGC
TGGGCTGAGCTTGGCCGACTGAACGGGTGCGGAAGTCTCAAGCCCTGAGTTCTCCAAG
ACTCTACAGCTTGAGTACCGCTTTGAGACAGTCCAGAAGCTACGCTTTGGAATCTATGAC
ATAGACAACAAGACGCCAGAGCTGAGGGATGATGACTTCTAGGGGGTGTGAGTGTCC
CTAGGACAGATTGTGCCAGCCAGTACTGACTCTCCCCTTGATGCTGAAGCCTGGAAAA
CCTGCTGGGCGGGGGACCATCACGGTCTCAGCTCAGGAATTAAGGACAATCGTGTAGTA
ACCATGGAGGTAGAGGCCAGAACTAGATAAGAAGGACTTCTGNGAAAATCAGATCCA
TTTCTGGAGTTCTTCCGCCAGGGTGTGGGAAATGGCACCTGGTGTACAGATCTGAGGTC
ATCAAGAACAACCTGAACCTACATGGAAGCGTTTCTCAGTCCCCGTTACGATTTCTGT
GGTGGGAACCCAGCACACCCATCCCAGTGCATGCTCCGATTATGACAGTGACGGNTCA
CATGATCTCATCGGTACCTTCCACACCAGCTTNGCCAGCTGCNAGCAGTCCCGGCTGAGT
TTTGATGCATCCACCCTGAGAAGCGCAGAAAAGAAAGCTACAGACTCTGGACTATNCGGT
A
    
```

Restriction Sites:

NotI-NotI

ACCN:

NM_152927

Insert Size:

2050 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_152927.1 , NP_690904.1
RefSeq Size:	2107 bp
RefSeq ORF:	1614 bp
Locus ID:	8904
UniProt ID:	Q99829
Cytogenetics:	20q11.22
Gene Summary:	<p>Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene encodes a calcium-dependent protein that also contains two N-terminal type II C2 domains and an integrin A domain-like sequence in the C-terminus. However, the encoded protein does not contain a predicted signal sequence or transmembrane domains. This protein has a broad tissue distribution and it may function in membrane trafficking. This gene and the gene for RNA binding motif protein 12 overlap at map location 20q11.21. Alternate splicing results in multiple transcript variants encoding different proteins. [provided by RefSeq, Aug 2008]</p> <p>Transcript Variant: This variant (4) uses an additional exon in its 5' UTR, compared to variant 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>