

Product datasheet for SC110251

DCTN3 (NM_024348) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DCTN3 (NM_024348) Human Untagged Clone
Tag:	Tag Free
Symbol:	DCTN3
Synonyms:	DCTN-22; DCTN22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC110251 representing NM_024348. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCCGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGGTCTGACTGACTTGCAGCGGCTACAGGCCGAGTGGAAGAGCTGGAGCGCTGGGTGTACGGG
CCGGGCGGGGCGCGCGGCTCACGGAAGGTGGCTGACGGCCTGGTCAAGGTGCAGGTGGCTTTGGGGAAC
ATTTCCAGCAAGAGGGAGAGGGTGAAGATTCTCTACAAAAGATTGAAGATCTGATCAAGTACCTGGAT
CCTGAGTACATCGACCGCATTGCCATACCTGATGCCTCTAAGCTGCAATTCATCCTAGCAGAGGAGCAG
TTTATCCTTTCCAGGTTGCACTCCTGGAGCAGGTGAATGCCTTGGTCCCATGCTGGACAGTGCTCAC
ATCAAAGCCGTTCTGAGCATGCTGCCCGCTGCAGCGCTTGGCCAGATCCACATTCAGCAGCAGGCT
CCATGGGGAGTGGGAGTCCGTGATGAGGCAGGAAGTTTAGTGGAAGATGTGGGCTTTGCCAGTTCCTT
TCTGTGCTACACTTTGGCCCTACAGGACCAAGTGTGTGGAATCACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_024348
Insert Size:	531 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024348.3
RefSeq Size:	954 bp
RefSeq ORF:	531 bp
Locus ID:	11258
UniProt ID:	O75935
Cytogenetics:	9p13.3
MW:	19.5 kDa
Gene Summary:	<p>This gene encodes the smallest subunit of dynactin, a macromolecular complex consisting of 10 subunits ranging in size from 22 to 150 kD. Dynactin binds to both microtubules and cytoplasmic dynein. It is involved in a diverse array of cellular functions, including ER-to-Golgi transport, the centripetal movement of lysosomes and endosomes, spindle formation, cytokinesis, chromosome movement, nuclear positioning, and axonogenesis. This subunit, like most other dynactin subunits, exists only as a part of the dynactin complex. It is primarily an alpha-helical protein with very little coiled coil, and binds directly to the largest subunit (p150) of dynactin. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]</p> <p>Transcript Variant: This variant (2) uses an alternate splice site in the 3' coding region, which results in a translational frameshift, compared to variant 1. The encoded isoform (2) has a shorter and distinct C-terminus, compared to isoform 1.</p>