

## Product datasheet for **SC111029**

### **SURF5 (MED22) (NM\_133640) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	SURF5 (MED22) (NM_133640) Human Untagged Clone
Tag:	Tag Free
Symbol:	SURF5
Synonyms:	MED24; SRB6; surf-5; SURF5
Mammalian Cell Selection:	Neomycin
Vector:	<u><a href="#">PCMV6-Neo</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)
Restriction Sites:	NotI-NotI
ACCN:	NM_133640
Insert Size:	1300 bp
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
--------------------	--



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_133640.3](#), [NP\\_598395.1](#)

**RefSeq Size:** 1448 bp

**RefSeq ORF:** 603 bp

**Locus ID:** 6837

**UniProt ID:** [Q15528](#)

**Cytogenetics:** 9q34.2

**Gene Summary:** This gene encodes a protein component of the mediator complex, which functions in the regulation of transcription by bridging interactions between gene-specific regulatory factors, RNA polymerase II, and general transcription factors. Alternatively spliced transcript variants encoding different isoforms have been observed. [provided by RefSeq, Jul 2013]  
Transcript Variant: This variant (b) represents the shorter transcript and encodes the longer isoform (b).