

Product datasheet for **SC111649**

LAIR2 (NM_002288) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LAIR2 (NM_002288) Human Untagged Clone
Tag:	Tag Free
Symbol:	LAIR2
Synonyms:	CD306
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_002288, the custom clone sequence may differ by one or more nucleotides

```
ATGTCTCCACACCTCACTGCTCTCCTGGGCCTAGTGCTCTGCCTGGCCCAGACCATCCACACGCAGGAGG
GGGCCCTTCCCAGACCTCCATCTCGGCTGAGCCAGGCACTGTGATCTCCCCGGGGAGCCATGTGACTTT
CATGTGCCGGGGCCCGTTGGGGTTCAAACATTCCGCCTGGAGAGGGAGGATAGAGCCAAGTACAAAGAT
AGTTATAATGTGTTTCGACTTGGTCCATCTGAGTCAGAGGCCAGATTCCACATTGACTCAGTAAGTGAAG
GAAATGCCGGGCTTTATCGCTGCCTCTATTATAAGCCCCCTGGATGGTCTGAGCACAGTGACTTCCTGGA
GCTGCTGGTGAAGAAAGCTCTGGAGGCCCGGACTCCCCGGACACAGAGCCCGGCTCCTCAGCTGGGACT
GTGCCAGGCACTGAAGCCTCCGATTTGATGCACCATGA
```

Restriction Sites:	NotI-NotI
ACCN:	NM_002288
Insert Size:	2770 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_002288.3](#), [NP_002279.2](#)

RefSeq Size: 702 bp

RefSeq ORF: 702 bp

Locus ID: 3904

UniProt ID: [Q6ISS4](#)

Cytogenetics: 19q13.42

Domains: IG

Protein Families: Secreted Protein

Gene Summary: The protein encoded by this gene is a member of the immunoglobulin superfamily. It was identified by its similarity to leukocyte-associated immunoglobulin-like receptor 1, a membrane-bound receptor that modulates innate immune response. The protein encoded by this locus is a soluble receptor that may play roles in both inhibition of collagen-induced platelet aggregation and vessel formation during placental implantation. This gene maps to a region of 19q13.4, termed the leukocyte receptor cluster, which contains 29 genes in the immunoglobulin superfamily. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Sep 2013]

Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).