

Product datasheet for **SC112397**

CCDC71 (NM_022903) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC71 (NM_022903) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC71
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_022903, the custom clone sequence may differ by one or more nucleotides

ATGAGTGTGGTGGTTCAGCATGTGGAGGAAAAAGCTGTGCACTCCTGGTCGCGCATCTCCACGGCAGGGA
AGAAGGCCCTGGAAGAGGCACTGTTGTCTTTAACCCAATGAGCCAGGATCTCAGTGCCACAGAGGCCCA
GCTTGTGGCCTTCCTGCAGGGCTGCGAGATGATGGCTTCCAACCTACCATCCTGCGCAGTGGTGATGTC
TATGGCTATAGTTCATGCACAGCTAATCCCCAAGCCAGACGAACTGCAAGCTCGTGCCCTAACCCAA
CTGCCACATCACCTCCAGCCAGTGCTCCCCGAAGTCCCATGCGGTTGCCTGCAGGTGGGCCACACTGCT
TCCCATGCCGCTATCTGGCAGACTGGCCAAAGCATCCACACCAGCCCTTGCCAAGCATGCTACCACCAAC
CTGCTGCTGAGCTCTCTGAAGCAATCAAGTGCCAGCCATGCCCGGGTGCAGCAGTGGGCTTCCCCACCC
ACCTTTATCCAGGTGTCTACCCTGCCATGCGGCTCTCTGTTGTCTTGAGGCCCTGGTTCCACTTAAGAC
TCCCATGCCCTGCTTGGGTGCCAAGCACAGGCACAGTCACTGCAGCTCTCACTTGCAGACTCTCCTCTG
AAACTGCGGAAAAGTTCAGGGAAGGGTCCGGGGAACCTCGGCCCAAAGCTCCAGAAAAACCACAAGCA
AGGGCCCCAAGTGTCTGACTCGCAAAGGCCCTGGGGCTGGACCCGACGAGGCTCTGGGCACCAGAGCAA
AACCAACAGAGCCACTGGGTCCCCAGTGTCCGGCGAATGAAAGGGGGCTCTGCCCTGGGCACCAAAACA
GCCCAGGCCAAGGTAGCTCGAAGCTGGCCAAAGCTGCTCGTCCCAAGCCAAGGTGGCTCGAACACAGG
CCAAGGTGCTAAGGCCGGGCCAAGGCCAAGGCAGCACAGGTCAAGGCTAAGGCCAAAGCTAAGGCAGC
ACAGGTCAAGGCCAAGGCCAAAGTCATGGCAGCATGGGCCAAGGCCAAGGCTAAAGCCAAGGCAGTACGG
GCCAAGGCCAAGGTGGCTCGGACCCAGCCAGGGGCAGAGGCAGGCCAAAGGGATCTGCTAAGGCCAGAA
CTACAAGGAAGGGCCAGAAAAACCGCCCTGAGACTGTTGGGCAGAAGAGGAAAAGGGCTGAGGAGGCAAA
AGATCTTCCTCCCAAGAAGAGAACACGGCTTGGCCCCGATCTCCTAAGGCATGGCTAGGGCCTGGAACA
GCAAAGCTGCTGAAGTTCCTGTCATAAAGGTAGATAGGCGGTCTCGGATGATGAGGTGCGGCAGCGGG
CTCAGCGGATCTCCGCTGAACCTGTACCTGTAATACGGCTCCAGCCATTGCTGCCATATTACAGCAGT
CTGA

Restriction Sites: ECoRI-NOT



[View online »](#)

ACCN:	NM_022903
Insert Size:	1380 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_022903.2</u> , <u>NP_075054.2</u>
RefSeq Size:	1820 bp
RefSeq ORF:	1404 bp
Locus ID:	64925
UniProt ID:	<u>Q8IV32</u>
Cytogenetics:	3p21.31