

Product datasheet for **SC112770**

PGAP3 (NM_033419) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PGAP3 (NM_033419) Human Untagged Clone
Tag:	Tag Free
Symbol:	PGAP3
Synonyms:	AGLA546; CAB2; hCOS16; PERLD1; PP1498
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_033419, the custom clone sequence may differ by one or more nucleotides

```
ATGGCCGGCCTGGCGGCGCGGTTGGTCCTGCTAGCTGGGGCAGCGGCGCTGGCGAGCGGCTCCCAGGGCG  
ACCGTGAGCCGGTGTACCGGACTGCGTACTGCACTGCGAAGAGCAGAAGTCTCTGGGGGCGCTCTGAA  
TCACTTCCGCTCCCGCCAGCCAATCTACATGAGTCTAGCAGGCTGGACCTGTCGGGACGACTGTAAGTAT  
GAGTGTATGTGGGTACCGTTGGGCTCTACCTCCAGGAAGGTCAAAAAGTGCCTCAGTTCATGGCAAGT  
GGCCCTTCTCCCGTTCTGTTCTTTCAAGAGCCGGCATCGGCCGTGGCCTCGTTTCTCAATGGCCTGGC  
CAGCCTGGTGATGCTCTGCCGCTACCGCACCTTCGTGCCAGCCTCTCCCCATGTACCACCTGTGTG  
GCCTTCGCTGGGTGCCCTCAATGCATGTTCTGGTCCACAGTTTTCCACACCAGGGACTGACCTCA  
CAGAGAAAATGGACTACTTCTGTGCCTCCACTGTCATCCTACACTCAATCTACCTGTGCTGCGTCAGGAC  
CGTGGGGCTGCAGACCCAGCTGTGGTCACTGCCTTCCGGGCTCTCCTGTGCTCATGCTGACCGTGCAC  
GTCTCCTACCTGAGCCTCATCCGCTTCGACTATGGCTACAACCTGGTGGCCAACGTGGCTATTGGCCTGG  
TCAACGTGGTGTGGTGGCTGGCCTGGTGCCTGTGGAACCAAGCGGGCGGCTGCCTCACGTGCGCAAGTGCCT  
GGTGGTGGTCTTGCTGCTGCAGGGGCTGTCCCTGCTCGAGCTGCTTGACTCCCACCGCTCTTCTGGGTC  
CTGGATGCCATGCCATCTGGCACATCAGCACCATCCCTGTCCACGTCTCTTTTTCAGCTTTCTGGAAG  
ATGACAGCCTGTACCTGCTGAAGGAATCAGAGGACAAGTCAAGCTGGACTGA
```



[View online »](#)

5' Read Nucleotide Sequence:	<p>>OriGene 5' read for NM_033419 unedited</p> <pre>CTATAGGGCGCCGCAATTCGGCACGAGGGGCGGCGAGCCAGGGAGAAAGGATGGCCGG CCTGGCGGCGGGTTGGTCTGTAGCTGGGGCAGCGGCGCTGGCGAGCGGCTCCAGGG CGACCGTGAGCCGGTGTACCGCGACTGCGTACTGCAGTGCGAAGAGCAGAACTGCTCTGG GGGCGCTCTGAATCACTCCGCTCCCGCCAGCCAATCTACATGAGTCTAGCAGGCTGGAC CTGTCCGGACGACTGTAAGTATGAGTGTATGTGGGTACCGTTGGGCTCTACCTCCAGGA AGGTCACAAAGTGCCTCAGTCCATGGCAAGTGGCCCTTCTCCCGTTCCCTGTTCTTTCA AGAGCCGGCATCGGCCGTGGCCTCGTTTTCTCAATGGCCTGGCCAGCCTGGTGATGCTCTG CCGCTACCGCACCTTCGTGCCAGCCTCCTCCCCAAGTTACCACACCTGTGTGGCCTTCGC CTGGGTGTCCCTCAATGCATGGTCTGGTCCACAGTTTTCCACACCAGGGACACTGACCT CACAGAGAAAATGGACTACTTCTGTGCCTCCACTGTCATCCTACACTCAATCTACCTGTG CTGCGTCAGGACCGTGGGGCTGCAGCACCCAGCTGTGGTCAAGTGCCTTCCGGGCTCTNCT GCTGCTCATGCTGACCGTGCACGCTCCTACCTGAGCCTCATNCGCTTCGACTATGGCTA CAACCTGGTGGGCCACGTGGCTATTGGGCTGGTCAACGTGGNTGTGGTGGCTGGCCTGG</pre>
3' Read Nucleotide Sequence:	<p>>OriGene 3' read for NM_033419 unedited</p> <pre>GCTATGGACCGGCGCCGCAATCTANGATCGAGTTTTTTTTTTTTTTTTTTTTTTTGGAGCTTGG AATTGATTTATTAAGCACATCCCTTGCCACCCTCCACCTTAAAAGTTTTTCAGTATCAAA AGAAGCTAGCGCAGGCCACCCGAGTCCCTGAGGTTGGAGTCTAGCATAGCTCCCTCCC TCAAAGAGGGACAAGGGGTGAGGGGAGAGCAAAAATCCAGTCTGCTTCAACCACGGAGA CTGCCTTTGGGATGGAAGTTTTCTGGAGTCCCTCCATTCTATTCTGTGGGGCAGGAAC ATGCCAGGGCTGTGGTAAATGGCAGGGTCCCTTTACCAGGGCGCAGGCATAGTGTGG CCCCTGCCTGCCCTGGGGCCACCCTGGAACAGTAAGACCCTCTGGGAAGGGAAGAAACCA GAGCCAGCCATAAGCTTTCTCCCTTCTCCCTCAAATCGCCCTTCTCTCTCCCAAAAT CAGGCTATGGCAAAGGGGCTGGCAGCAGGGGAATGGAGTTTTTGGAGGGGAAGGCCAGCC ATGCGAAGGCAAGGAGAAAAGGAGGATGGAGAGAGCCCAAGCCATCCCCATTTCTCTTCC TCTGCCTTCGAGGGGCTGCTTTTGGAGGGCACTGTTGTGTGGGGACCCTGTGATGTTTGC AGGAGGCGGGCCAAAGGGGACCCAGACAGGGGACATGCACACAGTCAGTCTGCCCTGG AGAGCTCCTTGGCTGTCCGCCAGGTGGCCGCCCTTGGCAAAGCTGGCGCTCCACAAGGG ATTATTGGTGGGAAAGGATTCCCTTTACCTGTTTCAAACCTTAACTGGGGCCCATTTCC CCGGATGTTCTCATTAAAATGCCACTAAACCGCCACCCACTTCCCTTACCAGCTTA ACCCCCACCCCTGCAACACCACTTTTGAGACCCCTTTTCCCGAAACCACCGTTTCCC ACTGGCTAT</pre>
Restriction Sites:	ECORI-NOT
ACCN:	NM_033419
Insert Size:	2620 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_033419.3](#), [NP_219487.3](#)

RefSeq Size: 2708 bp

RefSeq ORF: 963 bp

Locus ID: 93210

UniProt ID: [Q96FM1](#)

Cytogenetics: 17q12

Domains: Per1

Protein Families: Transmembrane

Gene Summary: This gene encodes a glycosylphosphatidylinositol (GPI)-specific phospholipase that primarily localizes to the Golgi apparatus. This ubiquitously expressed gene is predicted to encode a seven-transmembrane protein that removes unsaturated fatty acids from the sn-2 position of GPI. The remodeling of the constituent fatty acids on GPI is thought to be important for the proper association between GPI-anchored proteins and lipid rafts. The tethering of proteins to plasma membranes via posttranslational GPI-anchoring is thought to play a role in protein sorting and trafficking. Mutations in this gene cause an autosomal recessive form of neurologic hyperphosphatasia with cognitive disability (HPMRS4). Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2017]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).