

Product datasheet for **SC112796**

MRPS14 (NM_022100) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPS14 (NM_022100) Human Untagged Clone
Tag:	Tag Free
Symbol:	MRPS14
Synonyms:	COXPD38; DJ262D12.2; HSMRPS14; MRP-S14; S14mt
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_022100, the custom clone sequence may differ by one or more nucleotides

ATGGCGCCTTCATGCTGGGCTCGCTGCTGCGGACGTTCAAGCAGATGGTTCCTTCATCAGCTTCAGGCC
AAGTTCGAAGTCACTATGTAGACTGGAGAATGTGGCGGATGTGAAGAGACGAAAAATGGCCTATGAATA
CGCAGATGAGAGGCTACGTATTAATCACTCAGGAAGAATACCATTTGCCAAAAATCTTCAGGATGTG
GCTGATGAAGAAATTGCTGCCCTCCCCGGGATAGCTGTCCTGTTAGAATCAGAAATCGGTGTGTTATGA
CGTCCCGTCCGCGTGGTGTGAAGCGGCGCTGGAGGCTTAGTCGTATAGTCTCCGTCACCTAGCTGACCA
TGGGCAACTTTCTGGGATCCAGCGAGCGACATGGTAA



[View online »](#)

5' Read Nucleotide Sequence:

>OriGene 5' read for NM_022100 unedited
 AGCAATTGTATACGACTCATATAGGCGGCCGCGAATTCGCACGAGCCTCGTGCCGAATT
 CGGCACGAGGGGACAACATGGCGGCTTCATGCTGGGCTCGCTGCTGCGGACGTTCAAGC
 AGGTGAGGCTCTACTTTTATCCACACCGCTACCCCTACCGCTGAATCTCATAACCCAC
 GGGTCCCCACGGCTCTGCCACCTGTAGTGCTCGATGGTTCCTTCATCAGCTTCAGGCCA
 AGTTCGAAGTCACTATGTAGACTGGAGAATGTGGCGCGATGTGAAGAGACGAAAAATGGC
 CTATGAATACGCAGATGAGAGGCTACGTATTAATCACTCAGGAAGAATACCATTTTGCC
 AAAAATTCCTTCAGGATGTGGCTGATGAAGAAATTGCTGCCCTCCCCGGGATAGCTGTCC
 TGTTAGAATCAGAAATCGGTGTGTATGACGTCCCGTCCGCGTGGTGTGAAGCGCGCTG
 GAGGCTTAGTCGTATAGTCTTCCGTCACTAGCTGACCATGGGCAACTTTCTGGGATCCA
 GCGAGCGACATGGTAAATGAGCTCCAGAACCTATTGAGCTTGACGGGAAGCCAAGCTTGC
 AGTTCAGCAAGCAAAGATTTTTTTAATAGACCAAACCTAATCTCTACAGGGGCCAG
 TACAGTTGTTGGCTACCTGATGCTATCTCTAACTACTTTTAAATGAAGACATTTGG
 GTGTTTCATGTCAGTGAATTATCTTTCTTGGACTTAACANACATACTGNTCCCATAC
 ATAGTACANAGAGTAGCATAATAAATATTTCTGTNGAAAAAAAAAAAAAACTCGACTC
 TAGATTGCGGGCGCGGTATAGCTGGTTCCTGACAGATCCCGGTGGCATCCCTGTGACC
 CCTNCCCATGCCTCTCCTGGCCNTGNNAGTGCCACTCCATGCC

3' Read Nucleotide Sequence:

>OriGene 3' read for NM_022100 unedited
 CCGCGGGCCGAATCTANAATCGAGTTTTTTTTTTTTTTTTTTCACAGAAAAATTTATTA
 TGCTACTCTTTGTACTATGTATGGGAACAGTATGTTTGTTAAGTCCAAAAAAAAGATAA
 TTCCTGACATGAAACACCCAAATGTCTTCATTTTAAAGTAGTTTAAAGATAGCATCAG
 GTAGGCCAAACAAGTGTACTGGGCCCTGTAGAGATTAGGTTTGGTCTATTAATAAAAA
 TCTTTGCTTGCTGGAAGTCAAGCTTGGCTTCCTGCAAGCTCAATAGGTTCTGGAGCTC
 ATTTACCATGTCGCTCGCTGGATCCCAGAAAGTTGCCCATGGTCAGCTAAGTGACGGAAG
 ACTATACGACTAAGCCTCCAGCGCCGCTTCACACCACGCGGACGGGACGTCATAACACAC
 CGATTTCTGATTCTAACAGGACAGCTATCCCGGGGAGGGCAGCAATTTCTTCATCAGCC
 ACATCCTGAAGATTTTTGGCAAAATGGTATTCTTCTGAGTGAATTAATACGTAGCCTC
 TCATCTGCGTATTCATAGGCCATTTTTCGTCTCTTACATCGCGCCACATTCTCCAGTCTA
 CATAGTGACTTCGAACCTTGCCTGAAGCTGATGAAGGAACCATCGAGCACTACAGGTGGC
 AGAGCCGTGGGGACCCGTGGTTATGAGATTCAGCGGTGAGGGGTAGCCGGTGTGGATA
 AAAATAAAGCCCTGACCTGCTTGAACGTCCCGAGCAGCAGCCCAACATGAAGCCCGCATG
 TTGTCCCTTGTCGCGATTTCGCACGAGGTTGTCGCGAATTCCGGCCGCCCTAAGTGAGGC
 CTATTACAAATCTGACAGTTCATTAACCACTCTGTTATTTAACTCCCGTAACCCCTTCC
 GCCCATATGGTAACGGGCGGGTTTTACAACCTTTTGAAGAGC

Restriction Sites:

NotI-NotI

ACCN:

NM_022100

Insert Size:

820 bp

OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_022100.1, NP_071383.1</u>
RefSeq Size:	1081 bp
RefSeq ORF:	387 bp
Locus ID:	63931
UniProt ID:	<u>O60783</u>
Cytogenetics:	1q25.1
Domains:	Ribosomal_S14
Gene Summary:	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that belongs to the ribosomal protein S14P family. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2010] Transcript Variant: This variant (1) represents the shorter transcript and encodes the functional protein.