

Product datasheet for **SC113248**

HOXA5 (NM_019102) Human Untagged Clone

Product data:

| | |
|---------------------------|---------------------------------------------------------------------------------------------------|
| Product Type: | Expression Plasmids |
| Product Name: | HOXA5 (NM_019102) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | HOXA5 |
| Synonyms: | HOX1; HOX1.3; HOX1C |
| Mammalian Cell Selection: | None |
| Vector: | <u>pCMV6-XL4</u> |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| Fully Sequenced ORF: | >NCBI ORF sequence for NM_019102, the custom clone sequence may differ by one or more nucleotides |

```
ATGAGCTCTTATTTTGTAAACTCATTTTGCGGTCGCTATCCAAATGGCCCGGACTACCAAGTTGCATAATT
ATGGAGATCATAGTTCGGTGAGCGAGCAATTCAGGGACTCGGCGAGCATGCACTCCGGCAGGTACGGCTA
CGGCTACAATGGCATGGATCTCAGCGTCGGCCGCTCGGGCTCCGGCCACTTTGGCTCCGGAGAGCGCGCC
CGCAGCTACGCTGCCAGCGCCAGCGCGGCCCGCCGAGCCAGGTACAGCCAGCCGGCCACGTCCACGC
ACTCTCCTCAGCCCGATCCGCTGCCCTGCTCCGCCGTGGCCCCCTCGCCCGCAGCGACAGCCACCACGG
CGGGAAAACTCCCTAAGCAACTCCAGCGGCGCCTCGGCCGACGCCGCGCAGCACCCACATCAGCAGCAGA
GAGGGGTTGGCAGGCGTCCGGAGCCGAGGAGGACGCCCTGCCAGCAGCGAGCGAGCGAGTGGCGAGA
GCGAGCCGAGCCGGCGCGCCGCCCAACCCAGATCTACCCCTGGATGCGCAAGCTGCACATAAGTCA
TGACAACATAGGCGGCCGGAAGGCAAAAAGGGCCCGGACGGCCTACACGCGCTACCAGACCTGGAGCTG
GAGAAGGAGTTCCACTTCAACCGTTACCTGACCCGAGAAGGAGGATTGAAATAGCACATGCTCTTTGCC
TCTCCGAGAGACAAATTAATCTGTTTCCAAAACCGGAGAATGAAGTGGAAAAAAGATAATAAGCTGAA
AAGCATGAGCATGGCCGCGCAGGAGGGCCTTCGTCCTGA
```



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019102.2](#), [NP_061975.2](#)

RefSeq Size: 1332 bp

RefSeq ORF: 813 bp

Locus ID: 3202

UniProt ID: [P20719](#)

Cytogenetics: 7p15.2

Domains: homeobox

Protein Families: Transcription Factors

Gene Summary: In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. Methylation of this gene may result in the loss of its expression and, since the encoded protein upregulates the tumor suppressor p53, this protein may play an important role in tumorigenesis. [provided by RefSeq, Jul 2008]