

Product datasheet for **SC114101**

Apelin (APLN) (NM_017413) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Apelin (APLN) (NM_017413) Human Untagged Clone
Tag:	Tag Free
Symbol:	Apelin
Synonyms:	APEL; XNPEP2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL6</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_017413, the custom clone sequence may differ by one or more nucleotides

ATGAATCTGCGGCTCTGCGTGCAGGCGCTCCTGCTGCTCTGGCTCTCCTTGACCGCGGTGTGTGGAGGGT
CCCTGATGCCGCTTCCCGATGGGAATGGGCTGGAAGACGGCAATGTCCGCCACCTGGTGCAGCCCAGAGG
GTCAAGGAATGGGCCAGGGCCCTGGCAGGGAGGTTCGGAGGAAATCCGCCCCAGCGGCCCGCCTCTCC
CATAAGGGACCCATGCCTTTCTGA

5' Read Nucleotide Sequence:

>OriGene 5' read for NM_017413 unedited
CCCCCATACCCCGCCCGANGCCGCTTTGGGCGGTAGGCGTGTACGGTGGGAGGTCTAT
ATAAGCAGATCTCATTTAGGTGACACTATAGAATACAAGCTACTTGTCTTTTTGACGCG
GCCGCAATTCGGCACGAGGGGAGTTGCAGCCGCGGCCGAGCAGCCAGCCGCTAAGAAA
GAGCTCGCCGCTGCCGCTCCCGAGCCGCGGAGGCCAGCTTCGCGGCGCTGCCCCGCGC
GGGAGAGGAGGCTGCAGAAGAGCGGAGGCCAGCGGGAGCGGGGGCTCAGCGCGCA
CACTCAGCGGCCGGGAGCCTCCCGAGCTCCGCGCCCGCACGCGCCAGCCGCGGCTCGCG
CCTTTCTTGGCCTCCGGGCGCCGACCTCTCTCCCCGCGCCGCTCGCGGGGCGCGC
GCGGCCAAGGAGCAGCATGAATCTGCGGCTCTGCGTGCAGGCGCTCCTGCTGCTGTC
TCTCCTTGACCGCGGTGTGTGGAGGGTCCCTGATGCCGCTTCCCGATGGGAATGGGCTGG
AAGACGGCAATGTCCGCCACCTGGTGCAGCCCAGAGGGTCAAGGAATGGGCCAGGGCCCT
GGCAGGGAGGTTCGGAGGAAATCCGCGCCAGCGGCCCGCCCTCTCCATAAGGGACCCA
TGCCTTTCTGAAGCAGGACTGAAGGGGCCCCCAAGTGCCACCCCTGCGGTTATGTCTC
CTCCATAGAGTGGTCTGCTTCTCTGGAGGCCTCACGTCCATTAGCTCTCACCTCGCACC
TGCTGTAGCCACCAGTGGGCCAGCTCTTCTCACCTGCCTGCTTCCCCAGTGGCGTGCTC
CTGGCATGTAGTTTGGNATGAATCCCCGCTCTCACAGAAATCCGTCCAGTCCATCATTC
TGCTCC



[View online »](#)

3' Read Nucleotide Sequence:	<p>>OriGene 3' read for NM_017413 unedited TTTTTTAGATACTCTATGNACC GCGCCGAATCTANNATCGAGTTTTTTTTTTTTTTTT TTTATTATCAAATGTATTTATTGCTGAAAAACATAACATTTTTCAAGAAAGGCAAAGTGGC TAAAAAGTCTGAAAGTGTTCAAAGTAGATAAAATAGCAGAAGACACCCACCAAGGATA CAAAACAATTTTTAGTAGCGATCCTGCATTTAATCAGTATGTTCTTAAATAAACTGCTTT AAAAAATTTCTCAAATGACACTGCCAAAAAATTAGGACACCCAAACAGATGCCAGAAA ACCTGTAAGTGGGCTGGATTTTATGTGACCTGGTCATTAAGCATAGGGATTCAATTTGTG AGCTCATGGCAAATCAGTTTGAGGCCAGTTGACCTAGAACCATTGGGATGCAGGCCACT TACCTCCCTGCACGTGCAATATGTGGGCATGGGGACACCCATGACCCCCAGGGGAAAGGA ATGCCCTCTGAACATGACCTCCAAGAGTAAGGGCGAACTGTCAGCTTTTAACTGTTTA TTATTAAGACATATTTACACAGAACAATCTTTACAAACATTGAACACAGGGGAAGGGGAC AATTTCTTAATGAACAGGGCCTATATATCTTTGTATAAAATAGTATAAGAATCATAAACC ACCCCTTTAAATAGGGCACCCCTTAGCCACCCACTACTTTTTTTGTTCCCTAATCT CCCAATTTTCTTAAACATCCCCCACTTTTCTCCCTTTCCCAACGGGGCTTGGGCTTTG GCTTGCACGTTAATGGCAAGGTCGTTGAGGCCGGGAGCACAAACGGAAGGGGAGGGGAA AAGAAAGGCCCAATCCCCTGGAGGGGGCCGCATGGGGCGAAAGAGAATTTTTTCGGGA</p>
Restriction Sites:	NotI-NotI
ACCN:	NM_017413
Insert Size:	3500 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_017413.3 , NP_059109.3
RefSeq Size:	3125 bp

RefSeq ORF: 234 bp

Locus ID: 8862

UniProt ID: [Q9ULZ1](#)

Cytogenetics: Xq26.1

Protein Families: Druggable Genome, Secreted Protein

Gene Summary: This gene encodes a peptide that functions as an endogenous ligand for the G-protein coupled apelin receptor. The encoded preproprotein is proteolytically processed into biologically active C-terminal peptide fragments. These peptide fragments activate different tissue specific signaling pathways that regulate diverse biological functions including fluid homeostasis, cardiovascular function and insulin secretion. This protein also functions as a coreceptor for the human immunodeficiency virus 1. [provided by RefSeq, Feb 2016]