

Product datasheet for **SC114222**

NT5DC3 (NM_016575) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NT5DC3 (NM_016575) Human Untagged Clone
Tag:	Tag Free
Symbol:	NT5DC3
Synonyms:	GNN, FLJ11266, TU12B1TY, TU12B1-TY
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_016575, the custom clone sequence may differ by one or more nucleotides

ATGAGCAACTTGTGAATCCAGATGCCATTTTCTCAAACAATGAAATGAGCCTGTCAGACATTGAAATCT
ATGGCTTCGATTATGATTACACCTTGGTGTATTTATTCAAAGCACCTCCACACGCTGATATTTAATGCTGC
ACGGGACCTTCTCATCAATGAACACCGGTATCCAGCAGAAATCAGGAAGTATGAGTATGACCCAAATTTT
GCAATTCGTGGACTTCATTATGATGTACAGCGGGCAGTATTAATGAAGATCGATGCTTTTCATTATATCC
AGCTGGGAACTGTCTACAGAGGCCTCAGTGTGTCCCTGATGAAGAAATCATTGAAATGTACGAGGGGTC
CCACGTGCCCTTGGAGCAGATGAGTGACTTTTACGAAAGAGCTCTCATGAAACACGATGAAGCAGTTC
ATGGACATCTTCTCCCTGCCGAGATGACCCTCCTGCTGCGTGAATGAATACTTCTCAAGAACAACA
TCGACTATGAGCCTGTGCATCTGTACAAAGATGTCAAGGATTCAATTCGAGACGTCCACATCAAAGGAAT
AATGTACAGAGCAATTGAAGCAGACATTGAAAAGTACATCTGCTATGCTGAGCAGACCCGCGCAGTGTG
GCCAAACTGGCTGATCATGGCAAGAAGATGTTTCTCATACCAATAGCCCCAGTAGCTTTGTGGACAAAG
GGATGAGTTATATCGTTGGGAAAGACTGGAGGGACCTGTTTCGATGTGGTCATTGTTTCAGGCTGAGAAGCC
AACTTCTTTAATGATAAGCGGAGGCCTTTCCGAAAGATGAATGAGAAAGGTGCTTACTCTGGGATAAA
ATCCATAAGTTGCAGAAAGGCCAGATATACAAGCAGGGTAATTTATGAATTTTGAAGCTTACTGGAT
GGAGAGGATCCAGAGTGTGTATTTTGGTGACCATATATACAGTGACCTGGCGGATTTGACCCTAAAGCA
TGGCTGGAGGACTGGTGCAATCATCCAGAGTTGAGATCTGAGCTCAAAATCATGAACACGGAGCAATAC
ATTCAAACCATGACCTGGCTGCAGACCTTGACTGGCTTATTGGAACAGATGCAGGTTTCAGAGATGCTG
AGTCACAGCTGGTTTTCAGGAGTGGAAGAAAGGAAGGAGATGCGAGAAATGACCAAGAGTTTCTT
CAATGCCAGTTTGAAGCCTGTTCCGCACAGACCAGAACCCAACTACTTCTAAGCGCCTGTCGCGC
TTCGCTGACATCTACATGGCGTCTCTGAGCTGCCTCCTGAACTATGACGTGAGCCACACTTTCTACCCCC
GGAGGACTCCACTGCAGCACGAATGCCCGCTGGTCAGAAAGGCCCCCACTTCGGAACCCCTCTCT
GCAGGAGGCCAGGCCAAGTAG



[View online »](#)

Restriction Sites:	NotI-NotI
ACCN:	NM_016575
Insert Size:	2600 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_016575.1</u> , <u>NP_057659.1</u>
RefSeq Size:	7028 bp
RefSeq ORF:	1422 bp
Locus ID:	51559
Cytogenetics:	12q23.3