

Product datasheet for **SC115794**

SNRPD1 (NM_006938) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SNRPD1 (NM_006938) Human Untagged Clone
Tag:	Tag Free
Symbol:	SNRPD1
Synonyms:	HsT2456; Sm-D1; SMD1; SNRPD
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF sequence for NM_006938 edited ATGAAGCTCGTGAGATTTTTGATGAAATTGAGTCATGAAACTGTAACCATTGAATTGAAG AACGGAACACAGGTCCATGGAACAATCACAGGTGTGGATGTCAGCATGAATACACATCTT AAAGCTGTGAAAAATGACCCTGAAGAACAGAGAACCTGTACAGCTGGAAACGCTGAGTATT CGAGGAAATAACATTCGGTATTTTATTCTACCAGACAGTTTACCTCTGGATACACTACTT GTGGATGTTGAACCTAAGGTGAAATCTAAGAAAAGGGAAGCTGTTGCAGGAAGAGGCAGA GGAAGAGGAAGAGGAAGAGGACGTGGCCGTGGCAGAGGAAGAGGGGGTCTAGGCGATAA
5' Read Nucleotide Sequence:	>OriGene 5' read for NM_006938 unedited TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGGCAGGCGCCATGCGTCTG CAGTGCTCCGCGCGCTCTTGACGTCCGGAGCCCTGGAGTAGGCGCTTCCGGCCATTCAT ACTGCAGTCGGTCAAGTGTTCGGTTGAAGGATTCTGTGTGCTGTCGGACCCAGAGGGTGAC GGCGCCGTAGGATGAAGCTCGTGAGATTTTTGATGAAATTGAGTCATGAACTGTAACC ATTGAATTGAAGAACGGAACACAGGTCCATGGAACAATCACAGGTGTGGATGTCAGCATG AATACACATCTTAAAGCTGTGAAAAATGACCCTGAAGAACAGAGAACCTGTACAGCTGGAA ACGCTGAGTATTCGAGGAAATAACATTCGGTATTTTATTCTACCAGACAGTTTACCTCTG GATACACTACTTGTGGATGTTGAACCTAAGGTGAAATCTAAGAAAAGGGAAGCTGTTGCA GGAAGAGGCAGAGGAAGAGGAAGAGGACGTGGCCGTGGCAGAGGAAGAGGGGGT CCTAGGCGATAATGTCTCTCAAGATTTCAAAGTCATATGAGATTTGGGATATTTTTGT CAGGTTGTGTTTGTATGTCAGTTTTTATACACATAAATGTGGGACAGAGCTGTCTAT TTAGTATATCAAAGCTTTAGTAGTTTCCCTCCACATTACAGCATTACCACAGTGAGAGCT AAGCATTTCTACTGGGCAGTTTCATTNNTAGTTGATCAGGTTTTAAGTTTTTTGACTAAA TTTTTCCTTTCTCTATGATGAAAAGGTTAAAATAAACCTTACACAAATTAACAAAAA AAAAAAAAAAAAAAGTTGATTTATATTGCGCCCGGGTCTAATTGTTTCTGACAAATCCGT TGCAATCCTGG



[View online »](#)

3' Read Nucleotide Sequence: >OriGene 3' read for NM_006938 unedited
 CAATCTAGAGTCGAGTTTTTTTTTTTTTTTTTTTTTTTTTTTTTAATTAGGCTAAGGCTTTTA
 TTTTAACTTATTTCATCATAAAAAGAAAAAGAAAAATTTAGTTCAAAAACCTTAAACCT
 GATCAACTAAAAATGAAACTGCCAGTAAAAATGCTTAGCTCTCACTGTGGTAATTCGT
 GAATGTGGAGGAACTACTAAAACCTTTGATATACTAAATAGACAGCTCTGTCCCACATTT
 ATGTTTATTA AAAACTGACATAAAACAAACACAACCTGTACAAAAAATATCCCAAATCTCA
 TATGACTTTGAAATCTTGAGAGACATTATCGCCTAGGACCCCTCTTCCCTGCCACGGC
 CAGTCTCTTCCCTCTTCTCTTCTCTGCTCTTCTGCAACAGCTTCCCTTTTCTTAT
 ATTTACCTTAGTTCAACATCCACAAGTAGTGTATCCAGAGGTA AACTGTCTGGTAGAA
 TAAAAACCGAATGTTATTTCTCGAATACTCAGCGTTTCCAGCTGTACAGGTTCTCTGT
 TCTTCAGGGTCATTTTACAGCTTTAAGATGTGTATTCATGCTGACATCCACACCTGTGA
 TTGTTCCATGGACCTGTGTTCCGTTCTTCAATTCAATGTTACAAGTTCATGACTCATTT
 TCATCAAAAATCTCAGAGCTTCATCCTAGCGGCGCCGTACCCTCTGGGTCGACAGCA
 CACAGATCCCTTACCGAACCTGACCGATTGCATTATGATTGCCCGTAGCGCTACTCC
 AGGGGCTCCGNACGTAGTAACGCGCGGAGCACTGAATACGCATGGCGGNCCTCTGCCGA
 ATTTCCGCCCGCCTTAAGGGATCGAATAACAAAATCTGACGGTAACTTAAACGACCTTGC
 CTTTAGACCTNCCACGCCACCCCTACGCCCTTGGGTAACGGGCGGGGTTACCAATTTGG
 AAAGCCGTGGTTTTGGGCCAACAACTCCTTGACCCAAGCGGGA ACTTGA

Restriction Sites: NotI-NotI

ACCN: NM_006938

Insert Size: 780 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_006938.2](#), [NP_008869.1](#)

RefSeq Size: 1614 bp
RefSeq ORF: 360 bp
Locus ID: 6632

UniProt ID: [P62314](#)

Cytogenetics: 18q11.2

Domains: Sm

Protein Families: Druggable Genome, Stem cell - Pluripotency

Protein Pathways: Spliceosome, Systemic lupus erythematosus

Gene Summary: This gene encodes a small nuclear ribonucleoprotein that belongs to the SNRNP core protein family. The protein may act as a charged protein scaffold to promote SNRNP assembly or strengthen SNRNP-SNRNP interactions through nonspecific electrostatic contacts with RNA. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2014]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).