

Product datasheet for **SC117039**

CDK5 (NM_004935) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDK5 (NM_004935) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDK5
Synonyms:	LIS7; PSSALRE
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF sequence for NM_004935 edited ATGCAGAAATACGAGAACTGGAAAAGATTGGGGAAGGCACCTACGGAAGTGTGTTCAAG GCCAAAACCGGGAGACTCATGAGATCGTGGCTCTGAAACGGGTGAGGCTGGATGACGAT GATGAGGGTGTGCCGAGTTCCGCCCTCCGGGAGATCTGCCTACTCAAGGAGCTGAAGCAC AAGAACATCGTCAGGCTTCATGACGTCCTGCACAGCGACAAGAAGCTGACTTTGGTTTTT GAATTCTGTGACCAGGACCTGAAGAAGTATTTTGACAGTTGCAATGGTGACCTCGATCCT GAGATTGTAAAGTCATTCTCTCCAGCTACTAAAAGGGCTGGGATTCTGTCATAGCCGC AATGTGCTACACAGGGACCTGAAGCCCCAGAACCTGCTAATAACAGGAATGGGGAGCTG AAATTGGCTGATTTTGGCCTGGCTCGAGCCTTTGGGATCCCGTCCGCTGTTACTCAGCT GAGGTGGTCACACTGTGGTACCGCCACCGGATGTCCTCTTTGGGGCCAAGCTGTACTCC ACGTCCATCGACATGTGGTCAGCCGGCTGCATCTTTGCAGAGCTGGCCAATGCTGGGCGG CCTCTTTTTCCCGCAATGATGTGATGACAGTTGAAGAGGATCTTCCGACTGCTGGGG ACGCCCCAGGAGCAGTGGCCCTCTATGACCAAGCTGCCAGACTATAAGCCCTATCCG ATGTACCCGGCCACAACATCCCTGGTGAACGTCGTGCCAAAACCTCAATGCCACAGGGAGG GATCTGCTGCAGAACCTTCTGAAGTGTAAACCCTGTCCAGCGTATCTCAGCAGAAGAGGCC CTGCAGCACCCCTACTTCTCCGACTTCTGTCCGCCCTAG



[View online »](#)

5' Read Nucleotide Sequence:

>OriGene 5' read for NM_004935 unedited
 NGGTCAAATTTTGTATACGACTCCTATAGGGCGGCCGGAATTCGCACGAGGCGCAGGGG
 TCCCCGCGCCCGCCGAGCAGAATACGAGAACTGGAAAAGATGGGGAAGGCACCTACGG
 AACTGTGTTCAAGGCCAAAAACCGGGAGACTCATGAGATCGTGGCTCTGAAACGGGTGAG
 GCTGGATGACGATGATGAGGGTGTGCCGAGTCCGCCCTCCGGGAGATCTGCCTACTCAA
 GGAGCTGAAGCACAAGAACATCGTCAGGCTTCATGACGTCCTGCACAGCGACAAGAAGCT
 GACTTTGGTTTTTGAATTCTGTGACCAGGACCTGAAGAAGTATTTTGACAGTTGCAATGG
 TGACCTCGATCCTGAGATTGTAAAGTCATTCTCTTCCAGCTACTAAAAGGGCTGGGATT
 CTGTTCATAGCCGCAATGTGCTACACAGGGACCTGAAGCCCCAGAACCTGCTAATAAACAG
 GAATGGGGAGCTGAAATTGGCTGATTTTGGCTGGCTCGAGCCTTTGGGATTCCCGTCCG
 CTGTTACTCAGCTGAGGTGGTCACACTGTGGTACCGCCACCGGATGTCCTCTTTGGGGC
 CAAGCTGTACTCCACGTCCATCGACATGTGGTCAGCCGGCTGCATCTTTCAGAGCTGGC
 CAATGCTGGGCGGCTCTTTTCCCGCAATGATGTCGATGACCAGTTGAAGAGGATCTT
 CCGACTGCTGGGACGCCACCGAGGAGCAGTGGCCCTCATGACCAAGCTGCCAGACTA
 TAAGCCCTATCCGATGTACCCGGCCACACATCCCTGGTGAACGTCGTGCCAAACTCATG
 CCACAGGNAGGGATCTGTGCAGAACCTTCTGAGTGTAACCCTGC

3' Read Nucleotide Sequence:

>OriGene 3' read for NM_004935 unedited
 TATGGACCGCGCCGCAATCTAGGATCGAGTTTTTTTTTTTTTTTTTTTGTGGAAAGGA
 GCCAATTTATGAAATTAATAAAGTCCACAAAGGGGAGTGAGAAATTCGGNGCTCAGGG
 CACCCACACCCCGCGTGGGGCCAGCACCTGCTTGGAGCACAGACGCACCAGAGCACCC
 CCCACTTGTCTCACACCCTCTCAAGAGGGGGCCTTAAATAGGCCAGAGCCCAGCCCTGG
 AGGCCGGGGGTCCCGGGNCCTAGGGGCGGACACGAAGTCGNGAGAAGTAGGGGTGCGT
 GCAGGGCCCTCTTCTGCTGAGATACGCTGGACAGGGTTACACTTCAGAAGTTCTGCAGC
 AGATCCCTCCCTGTGGCATTGAGTTTGGGCACGACGTTACCAGGGATGTTGTGGCCGGG
 TACATCGGATAGGGCTTATAGTCTGGCAGCTTGGTCATAGAGGGCCACTGCTCCTCGGTG
 GGCGTCCCAGCAGTCGGAAGATCCTCTTCAACTGGTCATCGACATCATTGCCGGGAAAA
 AGAGGCCGCCAGCATTGGCCAGCTCTGCAAAGATGCAGCCGGCTGACCACATGTCGATG
 GACGTGGAGTACAGCTTGGCCCCAAGAGGACATCCGGTGGGCGGTACCACAGTGTGACC
 ACCTCAGCTGAGTAACAGCGGACGGGAATCCCAAGGCTCGAGCCAGGCCAAAATCAGCCA
 ATTTTCAGCTCCCATTCTGCTTATTAGCCGGTTCCTGGGCCTCAGGTCCTGTGTAGCC
 ACATGCGGCTATGACAGAATCCCACCCCTTTTAGTACCTGGAAGAGCATGACCTTACAAT
 CTCAGATCGAGGTCCCATTGCAACTGGCAAATACCTCTTCAGTCCCTGGTCCAGAATTCA
 AAACCAAGTCCCTTCTGCCCTGTGCAGGACTCCTGAAGCTGCCAC

Restriction Sites:

NotI-NotI

ACCN:

NM_004935

Insert Size:

1200 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004935.2](#), [NP_004926.1](#)

RefSeq Size: 1143 bp

RefSeq ORF: 879 bp

Locus ID: 1020

UniProt ID: [Q00535](#)

Cytogenetics: 7q36.1

Domains: pkinase, TyrKc, S_TKc

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Alzheimer's disease, Axon guidance

Gene Summary:

This gene encodes a proline-directed serine/threonine kinase that is a member of the cyclin-dependent kinase family of proteins. Unlike other members of the family, the protein encoded by this gene does not directly control cell cycle regulation. Instead the protein, which is predominantly expressed at high levels in mammalian postmitotic central nervous system neurons, functions in diverse processes such as synaptic plasticity and neuronal migration through phosphorylation of proteins required for cytoskeletal organization, endocytosis and exocytosis, and apoptosis. In humans, an allelic variant of the gene that results in undetectable levels of the protein has been associated with lethal autosomal recessive lissencephaly-7. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2015]

Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.