

Product datasheet for **SC118821**

Glutathione Peroxidase 4 (GPX4) (NM_002085) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Glutathione Peroxidase 4 (GPX4) (NM_002085) Human Untagged Clone
Symbol:	Glutathione Peroxidase 4
Synonyms:	GPx-4; GSHPx-4; MCSP; PHGPx; SMDS; snGPx; snPHGPx
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_002085, the custom clone sequence may differ by one or more nucleotides

```
ATGAGCCTCGGCCGCCTTTGCCGCTACTGAAGCCGGCGCTGCTCTGTGGGGCTCTGGCCGCGCCTGGCC  
TGGCCGGGACCATGTGCGCGTCCCGGGACGACTGGCGCTGTGCGCGCTCCATGCACGAGTTTTCCGCCAA  
GGACATCGACGGGCACATGGTTAACCTGGACAAGTACCGGGGCTTCGTGTGCATCGTCACCAACGTGGCC  
TCCCAGTGAGGCAAGACCGAAGTAACTACACTCAGCTCGTCGACCTGCACGCCGATACGCTGAGTGTG  
GTTTGGGATCCTGGCCTTCCCGTGTAAACAGTTCGGGAAGCAGGAGCCAGGGAGTAACGAAGAGATCAA  
AGAGTTCGCCCGGGCTACAACGTCAAATTCGATATGTTTCAGCAAGATCTGCGTGAACGGGACGACGCC  
CACCCGCTGTGGAAGTGGATGAAGATCCAACCAAGGGCAAGGGCATCCTGGGAAATGCCATCAAGTGGA  
ACTTCACCAAGTTCCTCATCGACAAGAACGGCTGCGTGGTGAAGCGCTACGGACCCATGGAGGAGCCCT  
GGTGATAGAGAAGGACCTGCCCCACTATTTCTAG
```



[View online »](#)

5' Read Nucleotide Sequence:	<p>>OriGene 5' read for NM_002085 unedited</p> <pre> CCCCCGATAGGATTTGTATACGACTTATATAGGCGGCCGCGCAATTCGCACGAGGGTGGG GGAGCCAGGAGGGGCGGGCGACGGGCGGGTATGGGCCGCGGGCGCAGGCTCCCCGGG CGCCGACAGCAGCGGTGCCAGAGCCGGGCGAGGCGGGCGCGGAGCCCTCGGCGGGCG AAGGCCCCAGCGTGCAGGCGCAGGAGGGCGGGCGCGCGGCGGAAGAAGCCCTGTCCCCG AGCTTGCACCCGAGATCCACGAATGTCCCAAGTCCCAGGACCCGTGCGCGTCCCGGGAC GACTGGCGCTGTGCGCGCTCCATGCACGAGTTTTCCGCCAAGGACATCGACGGGCACATG GTTAACCTGGACAAGTACCGGGGCTTCGTGTGCATCGTACCAACGTGGCCTCCCAGTGA GGCAAGACCGAAGTAAACTACACTCAGCTCGTGCACCTGCACGCCGATACGCTGAGTGT GGTTTGGGATCCTGGCCTTCCCGTGTAAACAGTTCGGGAAGCAGGAGCCAGGGAGTAAC GAAGAGATCAAAGAGTTCGCCGCGGGCTACAACGTCAAATTCGATATGTTACAGCAAGATC TGCGTGAACGGGACGACGCCACCCGCTGTGAAAGTGGATGAAGATCCAACCAAGGGC AAGGGCATCCTGTGAAATGCCATCAAGTGAACCTCACCAAGTTCCTCATCGACAAGACG GCTGCGTGGTGAAGCGCTACTGACCCATGGAGGAGCCCTGGTGTAGAGAAGGACCTGC CCCCTATTTCTAGCTCCACAAGTGTGTGGCCCCGCCGAGCCCTGCCACGCTCTTGG AGCCTTTCACCGCACTCATGACGGCCTGCCTGCCAACCTGCTGGTGGGGCAGACCCGA </pre>
3' Read Nucleotide Sequence:	<p>>OriGene 3' read for NM_002085 unedited</p> <pre> TTTATTCACCAAGTAGCCAGGGGTGGGGCGCCGAGCCAAGCCCAGCAGGCCATGGGAC CTTCTCCGGCGGGTGCACGCTGGATTTTCGGGTCTGCCCCACCAGCAGGTTTGCAGGC AGGCCGTATGAGTGCCGGTGAAGGCTCCAAGGGCGTGGGCAGGGGCTCGGGCGGGGCC ACACACTTGTGGAGCTAGAAATAGTGGGGCAGGTCCTTCTATCACCAGGGGCTCCTCC ATGGGTCCGTAGCGCTTACCACGCAGCCGTTCTTGTGATGAGGAACCTGGTGAAGTTC CACTTGATGGCATTTCACAGGATGCCCTTGCCCTTGGGTTGGATCTTCATCCACTCCAC AGCGGGTGGGCGTGTCCCGTTACGCAGATCTTGTGAACATATCGAATTTGACGTTG TAGCCCCGCGGCAACTCTTTGATCTTTCGTTACTCCCTGGCTCCTGCTTCCCGAACTGG TTACACGGGAAGGCCAGGATCCGCAAACCACTCAGCGTATCGGGCGTGCAGGTCGACG AGCTGAGTGTAGTTTACTTCGGTCTTGCCTCACTGGGAGGCCAGTGGTGACGATGCAC ACGAAGCCCCGGTACTTGTCCAGGTTAACCATGTGCCGTCGATGCTTGGCGGAAACT CGTGCATGGAGCGCGCACAGCGCCAGTCGTCTGGGACGCGCACGGNTCCTGNGACTTGG GACATCCGGGATCTCCGGTCCGAAGCTGCGGTGACAGGGCTTNTCCGCCGGCGCCGNG CCCTTCTGNNGTGCACGCTGGGGCCCTCCGCCGCCGAGGGGCTCGCGCGGNCNCTG CCCCGCTCTGGCACGNTGNCTGCNGCGCCGGGNNAGCTGCGCCGGGCGGCCATACCCG CCGTGTTCCGCCCTCTGGGTCCCCACCCTNCTGCGCATNCGGGCCGCTTATGGAGCGTT ATACAAATTTGACGGTCCACGAGCTTGTCTATTGACTCCCACGTTACCCTACCGNCA TTGCTG </pre>
Restriction Sites:	ECoRI-NOT
ACCN:	NM_002085
Insert Size:	1000 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). The expression of this clone is not guaranteed due to the nature of selenoproteins.
OTI Annotation:	This clone encodes a selenoprotein containing the rare amino acid selenocysteine (Sec). Sec is encoded by UGA codon, which normally signals translational termination. Expression of this clone is not guaranteed due to the nature of selenoproteins.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_002085.1</u> , <u>NP_002076.1</u>
RefSeq Size:	942 bp
Locus ID:	2879
UniProt ID:	<u>P36969</u>
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome
Protein Pathways:	Arachidonic acid metabolism, Glutathione metabolism
Gene Summary:	<p>The protein encoded by this gene belongs to the glutathione peroxidase family, members of which catalyze the reduction of hydrogen peroxide, organic hydroperoxides and lipid hydroperoxides, and thereby protect cells against oxidative damage. Several isozymes of this gene family exist in vertebrates, which vary in cellular location and substrate specificity. This isozyme has a high preference for lipid hydroperoxides and protects cells against membrane lipid peroxidation and cell death. It is also required for normal sperm development; thus, it has been identified as a 'moonlighting' protein because of its ability to serve dual functions as a peroxidase, as well as a structural protein in mature spermatozoa. Mutations in this gene are associated with Sedaghatian type of spondylometaphyseal dysplasia (SMDS). This isozyme is also a selenoprotein, containing the rare amino acid selenocysteine (Sec) at its active site. Sec is encoded by the UGA codon, which normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, designated the Sec insertion sequence (SECIS) element, that is necessary for the recognition of UGA as a Sec codon, rather than as a stop signal. Transcript variants resulting from alternative splicing or use of alternate promoters have been described to encode isoforms with different subcellular localization. [provided by RefSeq, Dec 2018]</p> <p>Transcript Variant: This variant (1) represents the predominant transcript and encodes the canonical isoform A. A similar isoform in rat (L-form) has been shown to be localized to the mitochondria (PMID:9988735). CCDS Note: This CCDS ID represents a selenoprotein that is functionally sensitive to selenium levels (e.g., PMID:14642406). Position 73 of the protein is a selenocysteine residue, which is cotranslationally incorporated when a UGA codon, which is normally a stop codon, is recognized by a specific selenocysteyl-tRNA in conjunction with a stem-loop SECIS element in the 3' UTR.</p>