

Product datasheet for SC120278

Selenoprotein M (SELENOM) (NM_080430) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Selenoprotein M (SELENOM) (NM_080430) Human Untagged Clone
Symbol:	Selenoprotein M
Synonyms:	SELM; SEPM
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL6</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF within SC120278 sequence for NM_080430 edited (data generated by NextGen Sequencing) ATGAGCCTCCTGTTGCCTCCGCTGGCGTGCTGCTGCTTCTCGCGGCGCTTGTGGCCCCA GCCACAGCCGCCACTGCCTACCGCCGGACTGGAACCGTCTGAGCGGCCTAACCCGCGCC CGGGTAGAGACCTGCGGGGGATGACAGCTGAACCGCCTAAAGGAGGTGAAGGCTTTCGTC ACGCAGGACATTCCATTCTATCACAACCTGGTGATGAAACACCTCCCTGGGGCCGACCCT GAGCTCGTGCTGCTGGGCCGCCGCTACGAGGAAGTAGAGCGCATCCCACTCAGTGAAATG ACCCGCAAGAGATCAATGCGCTAGTGACGAGAGCTCGGCTTCTACCGCAAGGCGGCGCCC GACGCGCAGGTGCCCCCGAGTACGTGTGGGCGCCCGCAAGCCCCAGAGGAACTTCG GACCACGCTGACCTGTAG Clone variation with respect to NM_080430.2 >OriGene 5' read for NM_080430 unedited GGCGGTAGGCGGTACGGTGGGAGGTCTATATAAGCAGAGCTCATTTAGGTGACACTATA GAATACAAGCTACTTGTCTTTTTGCAGCGGCCGCAATTCGGCACGAGGCTTGCGCATC ATCTGCGGCCGGGTCCCGATGAGCCTCCTGTTGCCTCCGCTGGCGCTGCTGCTTCTC GCGGCGCTTGTGGCCCCAGCCACAGCCGCACTGCCTACCGCCGGACTGGAACCGTCTG AGCGGCCTAACCCGCGCCCGGTAGAGACCTGCGGGGGATGACAGCTGAACCGCCTAAAG GAGGTGAAGGCTTTCGTACGCAGGACATTCCATTCTATCACAACCTGGTGATGAAACAC CTCCCTGGGGCCGACCCTGAGCTCGTGCTGCTGGGCCGCCGCTACGAGGAAGTAGAGCGC ATCCCACTCAGTGAAATGACCCGCAAGAGATCAATGCGTAGTGACGAGAGCTCGGCTTC TACCGCAAGGCGGCGCCGACGCGCAGGTGCCCCCGAGTACGTGTGGGCGCCCGCAAG CCCCCAGAGGAACTTCGACACGCTGACCTGTANGTCCGNGGCGGCGGAGCTGGG ACCTACCTGCCTGAGTCTGGAGACAGAATGAAGCGCTCAGCATCCCGGAATACTTCTC TTGCTGAGAGCCGATGCCCGTCCCCGGGCC
5' Read Nucleotide Sequence:	
Restriction Sites:	NotI-NotI
ACCN:	NM_080430


[View online »](#)

Insert Size:	800 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). The expression of this clone is not guaranteed due to the nature of selenoproteins.
OTI Annotation:	This clone encodes a selenoprotein containing the rare amino acid selenocysteine (Sec). Sec is encoded by UGA codon, which normally signals translational termination. Expression of this clone is not guaranteed due to the nature of selenoproteins.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_080430.2 , NP_536355.1
RefSeq Size:	739 bp
Locus ID:	140606
UniProt ID:	Q8WWX9
Cytogenetics:	22q12.2
Gene Summary:	The protein encoded by this gene belongs to the selenoprotein M/SEP15 family. The exact function of this protein is not known. It is localized in the perinuclear region, is highly expressed in the brain, and may be involved in neurodegenerative disorders. Transgenic mice with targeted deletion of this gene exhibit increased weight gain, suggesting a role for this gene in the regulation of body weight and energy metabolism. This protein is a selenoprotein, containing the rare amino acid selenocysteine (Sec). Sec is encoded by the UGA codon, which normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, designated the Sec insertion sequence (SECIS) element, that is necessary for the recognition of UGA as a Sec codon, rather than as a stop signal. [provided by RefSeq, Dec 2016]