

Product datasheet for **SC122366**

MRPL44 (BC012058) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPL44 (BC012058) Human Untagged Clone
Tag:	Tag Free
Symbol:	MRPL44
Synonyms:	39S ribosomal protein L44, mitochondrial; FLJ12701; FLJ13990; FLJ37688; L44MT; mitochondrial ribosomal protein L44; MRP-L44
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC012058 edited

```
GTTTTCTCTCGTGAATGGCGTCCGGCTGGTAAGATTGCTGCAGCAGGGACATCGCTGC
CTCCTGGTCCAGTCGCCCCAAGCTGGTCCCTCCGGTTCGGGGAGTGAAGAAGGGATTC
CGCGCCGCTTCCGCTTCCAGAAGGAGTTAGAGCGGCAGCGCCTTCTGCGGTGCCCGCC
CCGCCCCGTGCGCGTTCAGAGAAGCCGAAGTGGGATTACCATGCAGAAATACAAGCTTTT
GGACATCGGTTACAGGAAAACCTTTTCTTAGATCTTCTCAAAACTGCATTTGTTAATAGC
TGCTATATTAAGTGAAGGAGCCAAACGCCAACAACCTGGGATAGAGAAAGAAGCTGTT
CTTCTGAATCTTAAAGTAAATCAAGAACTATCCGAACAAGGGACATCTTTTTACAGACT
TGCCTTACACAGTTTCTTGAAGACGAGTACCCAGACATGCCCACTGAAGGCATAAAAAAT
CTTGTTGACTTTCTCACTGGTGAGGAAGTCGTGTGTACAGTGGCTAGAAACTGGCTGTG
GAGCAGTTAACTGAGTGAAGAATCCAGTCCCCAGCTGTGTTACAGCAGACTTTT
TTTGAGTTATTGGAGCCCTGTTACAGAGCAGTGGACCTGAGAGGACTGCACTTTTCATC
AGGGACTTCTTAATTACTCAAATGACTGGAAAAGAGCTCTTTGAGATGTGGAAGATAATA
AATCCCATGGGGCTATTGGTAGAAGAACTGAAGAAAAGGAATGTTTCAGCTCCTGAATCA
AGACTTACTAGGCACTGGTGGCACCACAGCTTTGCCTTTGTATTTTGTGGCTTATAC
TGTGATAAAAAGTTGATTGCAGAAGGACCTGGGAAAACAGTATTGGTTGCAGAAGAAGAG
GCTGCTCGAGTGGCCCTTAGAAAACCTTATGGATTACAGAAAATAGACGGCCGTGGAAC
TATTCCAAGCCCCAAGAAAACCTTGAGAGCAGAAAAGAGCATCACTGCCAGCTAGCCGCCA
TGGATGCAGCAGCCTGAACTTGAGAGCGAAAAGTGAATAAATGTCAAAGGTGTTTCAAG
CCAGACATTTTCACAATTGTGAAGAAAATAGATGTTTTGTTTCTGTTTTTACTGTGTTCC
CAAAATTAATAAATGTTAACCAAGTCAAAAAAAAAAAAAAAAAAAAA
```



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC012058 unedited AAAGAGCTGGGTTTTGTTTACGACTTCCTTGGGCGGCCGCAATTCGCACGAGGTTTTCTCTCGTGAATGGCGTCCGGGCTGGTAAGATTGCTGCAGCAGGGACATCGCTGCCTCCTGCTCCAGTCGCCCCAATCTGGTCCCTCCGGTTCGGGGAGTGAAGAAGGGATTCCGCGCCGCCTTCCGCTTCCAGAAGGAGTTAGAGCGGCAGCGCCTTCTGCGGTGCCCGCCGCCCGCCCGTTACAGGAAAACCTTTTCTTAGATCTTCTCAAACTGCATTTGTTAATAGCTGCTATATTAAGAGTGAGGAGCCAAACGCCAACAACTTGGGATAGAGAAAAGAAGCTGTTCTTCTGAATCTTAAAAGTAATCAAGAACTATCCGAACAAGGGACATCTTTTTTACAGACTTGCCTTACACAGTTTCTTGAAGACGAGTACCCAGACATGCCCACTGAAGGCATAAAAAATCTTGTGACTTTTCTACTGGTGAAGTACGTGTGTACGTGGCTAGAACTTGGCTGTGGAGCAGTTAACACTGAGTGAAGAATCCAGTGCCTCCAGCTGTGTACAGCAGACTTTCTTTGCA GTTATTGGAGCCCTGTTACAGAGCAGTGGACCTGAGAGGACTGCACCTTTCATCAGGGAC TTCTTAATTACTCANATGACTGGAAAAGAGCTCTTTGAGATGTGGAAGATAATAAATCCC ATGGGGGCTATTGGTAGAAGAAGTGAAGAANAAGGAATGTTTCAGCTCCTGATCAAGACT TACTANGCAGTCTGGTGGCACCACAGCTTGCCTTGTATTTTGTGGGCTATACTGAGAT AA
Restriction Sites:	Please inquire
ACCN:	BC012058
Insert Size:	1185 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	no
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC012058.1</u> , <u>AAH12058.1</u>
RefSeq Size:	1185 bp
Locus ID:	65080
Cytogenetics:	2q36.1

Gene Summary:

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. [provided by RefSeq, Jul 2008]