

Product datasheet for **SC122382**

TLE6 (BC007329) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TLE6 (BC007329) Human Untagged Clone
Tag:	Tag Free
Symbol:	TLE6
Synonyms:	GRG6; PREMBL
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC007329 edited

```
GGCAGGAGGCTTGGCTGACCTCCGCAGCTTCCGCCGACTCTGGCTAAAGTCTTGGAGGC
TACTGCCTTGAAGATGACCTCTAGGGACCAGCCAGACCCAAGGGCCCCGAAAAGCAC
TTTCGCTTGTCTGGGATCTCGAACTCTGAGAGCTCTCCGACGCTGAATTATCAGGGCAT
TCTAAATCGGCTCAAGCAGTTCCCCAGGTGAGAGCAATTCAGGGGCCGGAGGTCTCACG
CTGGGAAGGAGCCCCACCTGGCCTCTCTCCGGGGGTCTAGTGAGTGGTCTAGAGGCTG
CTGGCCCAAGAGCCAGGCTCAGGAACAGGAAGCCATGCAGATACCCAGCATGGAGAACC
GGCAATACCCATCCCAGCTTTTGATTAAACAAAATTACCAACAAATATTAATTTTGTAT
TTTCCTGTTCGAAAAAAAAAAAAAAAAAAAA
```

5' Read Nucleotide Sequence:

```
>OriGene 5' read for BC007329 unedited
TACGACTCACTATAGGGCGGCCGGAATTCGGCAGAGGCTTGGCTGACCTCCGCAGCTT
CCGCCGACTCTGGCTAAAGTCTTGGAGGCTACTGCCTTGAAGATGACCTCTAGGGACCA
GCCCAGACCCAAGGGCCCCCGAAAAGCACTTCGCCTTGTCTGGGATCTCGAACTCTGA
GAGCTCTCCGACGCTGAATTATCAGGGCATTCTAAATCGGCTCAAGCAGTTCCCCAGGTG
AGAGCAATTCAGGGGCCGGAGGTCTCACGCTGGGAAGGAGCCCCACCTGGCCTCTCTCC
CGGGGGTCTAGTGAGTGGTCTAGAGGCTGCTGGCCCAAGAGCCAGGCTCAGGAACAGGA
AGCCATGCAGATACCCAGCATGGAGAACC CGGCAATACCCATCCCAGCTTTTGATTAAAC
AAAATTACCAACAAATATTAATTTTGTATTTTCTGTTTCGAAAAAAAAAAAAAAAAAAC
TCGACTCTAGATTGCGGCCGCGGTCTAGCTGTTTCTGAACAGATCCCGGGTGGCATCC
CTGTGACCCCTCCCAGTGCCTCTCCTGGCCCTGGAAGTTGCCACTCCAGTGCCACCAG
CCTTGTCTAATAAAATTAAGTTGCATCATTTTGTCTGACTAGGTGTCCTTCTATAATAT
TATGGGTGGAAGGGGGTGGNNNNNNNNNGNNNNNGNNNNNNNNNNNNNNNNNNNNNN
GGNNNGGGGGNNNNNGGGNNNNNNNNNGGGGGGGGGNNNGNNNGGGGGNNNGG
GNGNNGGGNG
```

Restriction Sites: Please inquire



[View online »](#)

ACCN:	BC007329
Insert Size:	449 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC007329.1 , AAH07329.1
RefSeq Size:	449 bp
Locus ID:	84846
Cytogenetics:	19p13.3