

Product datasheet for **SC122591**

DIO3 (NM_001362) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DIO3 (NM_001362) Human Untagged Clone
Symbol:	DIO3
Synonyms:	5DIII; D3; DIOIII; TXDI3
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:

```
>OriGene sequence for NM_001362 edited
CGCCAGGCCACGTCGCGGTTGGTGGTCCGAGAGGGCGAGGGTCCCAGGGGGCTTCGGGG
CCTGCAGCCACCATGCTCCGCTCCCTGTGCTTCACTCCTTGAGGCTCTGCGCCAGACC
GCCTCGTGCCCTCGTGCTTCCCGCGCTTCTCGGCACGGCCTTCATGCTCTGGCTTCTC
GATTTCTGTGTATCCGCAAGCATTTCCTGGGCCGCCGCCCGGGGGCAGCCCGAGCCC
GAAGTGGAGCTCAACAGTGAAGGCGAGGAGGTGCCTCCCGATGACCCGCCCATCTGCGTG
TCCGACGACAACCGCCTGTGCACCCTGGCGTCGCTCAAGGCGGTGTGGCATGGCCAGAAG
TTGGATTTCTTCAAGCAGGCGCACGAGGGCGGTCCGGCGCCCAACTCCGAGGTGGTTCTG
CCCGACGGCTTCCAGAGCCAGCACATCCTCGACTACGCGCAAGGGAACCGCCCGTGGTT
CTCAATTCGCGAGCTGCACCTGACCACCGTTCATGGCGCGCATGAGCGCCTTCCAGCGC
CTGGTCACTAAGTACCAGCGCGAGCTGACTTCCCTCATCATCTACATCGAGGAAGCGCAC
CCCTCCGACGGCTGGGTACCACGGACTCTCCCTACATCATCCCACAGCACCGGAGCCTG
GAGGACCCGGTACGCGCAGCGAGGGTACTGCAGCAAGGTGCACCCGGCTGCGCTCTGGTC
CTCGACACCATGGCCAATCCAGCAGCTCGGCCTATGGCGCTACTTCGAGCGTCTCTAT
GTCATCCAGAGTGGCACTATTATGTACCAGGGCGCCGTGGCCCCGACGGCTACCAGGTC
TCTGAGCTGCGCACTTGGTTGGAACGCTATGATGAGCAACTGCACGGCGCTCGGCCCCGG
AGGGTGTAAACATCCAACGGACAATTGACTGAACTTGGTGGGCTGGGCCTTCGAGCCTTC
GAAGCCCACGTGCAAGCGCCTCAAACCAAGTACGCTTGGCGAGGCCCAAGTGACACTGA
TGTGCTGAGCCACCATTTAGACTGAGTCTGCACCCTCAGCCACATGAACAATCTCCCCT
ACCTCCCTGGGACTCTGCTTCTGTAAGTGTCTCATTACACCTGCCTGGCTCACTGGAAA
TCCTCTTTGAGCGCGGGATATGGCTTGGCCTTGTCCGTGTGCCCCAGGACTTTGCCTC
TACAGCATTTTCTTACACCCCTCCCGAGCGTCCCTCAGCCAAGTGCTTTGGCCCGGTG
CTTCCCGCAGCTGCACAGAGACCTTGGCCACGCCCGCGCCCTGAGCGCAGCTGGGTTT
CAGGAGACTCTCAGCTCAGCTAGCTAGTTGCCTGGCACCCACCTGTGCGCGCGGAGAG
GGGGTTCCCTGTTGCTTTTGTGCTGTTTCTGTCCCTGGTAGGGGAAGTATGTCGTGG
ATGGGGAGGGGTGGGCAGGGTAGTTTCCCGCTTGTGTTGGGTGCACAGGAGCCCACT
GCTGATGACGAATATCTCTAACTGGTCTTGACCACGAGCTAGTTCTGAATTGCAGGGGC
CTCAAAGCAGCACCTAACCTTGAGGGGGAGGGTGTCTGGGTTTCCGTGAGGTAACCAC
CTTAAATGGGAGGGAAGTTGGGGTGTCTGCTTTGGGACCAGAGGAAGATAGCTTGAGAGG
CATTGGCGAGGTTCCGACGCGCCAGGGAGAGAGAAAAAGCTGAGACTCCTGGGGATGA
CGTTGGGGTGATGGAGTCCGGGGAAAGAGAGGTGGGGGAGAGCCTGAGGTCCCCAAGTG
AGGGGAGTCTTAGCAGAGCTGCTGATTGTGGGCTGGGAGGTGGAGGGCCCTGATTCCG
AAGGCCATTTGGTGAAGTGTGCTGAAATATTTCTGTATATAAACTTCTTCAATCT
ACAATAATAAAGGCTTGAGGTAACTGCTTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```

5' Read Nucleotide Sequence:	>OriGene 5' read for NM_001362 unedited NGGGGATAGGATTTGTATACGACTTATATAGGCGGCCGCGCATTCCCGGGATATCGTCGA CCCACGCGTCCGCGCCAGGCCACGTGCGGGTTGGTGGTCGGAGAGGGCGAGGGGTCCCAG GGGGCTTCGGGGCCTGCAGCCACCATGCTCCGCTCCCTGCTGCTCACTCCTTGAGGCTC TGCGCCCAGACCGCTCGTGCCTCGTCTTCCCGCGCTTCTCGGCACGGCCTTCATG CTCTGGCTTCTCGATTTCTGTGTATCCGCAAGCATTTCCTGGGCCGCCGCCGGGGG CAGCCCGAGCCCGAAGTGGAGCTCAACAGTGAAGGCGAGGAGGTGCCTCCCGATGACCCG CCCATCTGCGTGTCCGACGACAACCGCCTGTGCACCCTGGCGTCGCTCAAGGCGGTGTGG CATGGCCAGAAGTTGGATTTCTCAAGCAGGCGCACGAGGGCGGTCCGGCGCCAACTCC GAGGTGGTTCTGCCGACGGCTTCCAGAGCCAGCACATCCTCGACTACGCGCAAGGGAAC CGCCCGTGGTTCTCAATTTCCGCGAGCTGCACCTGACCACCGTTCATGGCGCGCATGAGC GCCTTCCAGCGCCTGGTCACTAAGTACCAGCGCAGCTCGACTTCTCATCATCTACATC GAGGAAGCGCACCCCTCCGACGGCTGGGTACCACGGACTCTCCCTACATCATCCCACAG CACCGGAGCCTGGAGACCGGGTCAGCGCAGCGAGGGTACTGCAGCAAGGTGCACCCGGC TGGCTCTGGTCTCGACACCATGGCCAACTCCAGCAGCTCGGCCTATGGCGCCTACTTC GAGCGTCTCTATGTCATCCAGAGTGGCACTATTATGTACCCAGGCGGCCGTGGCCCGAC GGCTN
Restriction Sites:	Please inquire
ACCN:	NM_001362
Insert Size:	2098 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). The expression of this clone is not guaranteed due to the nature of selenoproteins.
OTI Annotation:	This clone encodes a selenoprotein containing the rare amino acid selenocysteine (Sec). Sec is encoded by UGA codon, which normally signals translational termination. Expression of this clone is not guaranteed due to the nature of selenoproteins.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001362.2</u> , <u>NP_001353.3</u>
RefSeq Size:	2120 bp
Locus ID:	1735
UniProt ID:	<u>P55073</u>
Cytogenetics:	14q32.31

Protein Families: Druggable Genome

Gene Summary: The protein encoded by this intronless gene belongs to the iodothyronine deiodinase family. It catalyzes the inactivation of thyroid hormone by inner ring deiodination of the prohormone thyroxine (T4) and the bioactive hormone 3,3',5-triiodothyronine (T3) to inactive metabolites, 3,3',5'-triiodothyronine (RT3) and 3,3'-diiodothyronine (T2), respectively. This enzyme is highly expressed in pregnant uterus, placenta, fetal and neonatal tissues, and thought to prevent premature exposure of developing fetal tissues to adult levels of thyroid hormones. It regulates circulating fetal thyroid hormone concentrations, and thus plays a critical role in mammalian development. Knockout mice lacking this gene exhibit abnormalities related to development and reproduction, and increased activity of this enzyme in infants with hemangiomas causes severe hypothyroidism. This protein is a selenoprotein, containing the rare selenocysteine (Sec) amino acid at its active site. Sec is encoded by the UGA codon, which normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, designated the Sec insertion sequence (SECIS) element, that is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. [provided by RefSeq, May 2016]