

Product datasheet for SC122947

C20orf195 (FNDC11) (NM_024059) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C20orf195 (FNDC11) (NM_024059) Human Untagged Clone
Tag:	Tag Free
Symbol:	C20orf195
Synonyms:	C20orf195
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_024059 edited

```

GGAAGAGGCCCCAGCACTGACCTCCGTGGGGGTGGAGATGAGGAGGATGGAAGGGTGTC
TTCCTCCAGCATCTTCTGAAGCTCCCGGATAATGAGCACCCATGTGGCAGGCCTGGGCC
TGGACAAGATGAAGCTGGGCAATCCCCAGTCCTTCTGGACCAGGAGGAGGCAGATGACC
AGCAGCTGCTGGAACCAGAGGCGTGGAAGACCTACACCGAGCGCCGAATGCCCTGCGTG
AGTTCTGACCTCGGACCTGAGTCCGCACCTGCTCAAGCGCCACCACGCCCGCATGCAGC
TGCTGCGTAAGTGCTCCTACTACATCGAGGTCCTGCCAAGCACCTGGCCCTGGGCGACC
AGAACCCGCTGGTGCTGCCTAGCGCCTTGTTCCAGCTCATCGACCCCTGGAAGTTCAGC
GCATGAAGAAGGTGGGCACAGCTCAGACCAAGATCCAGCTCCTGCTGCTCGGGGACCTGT
TGGAACAGCTCGACCATGGCCGTGCTGAGCTGGATGCCCTGCTCCGGTCGCCAGACCCAC
GGCCCTTCTGGCCGACTGGGCGCTGGTGGAGCGCGGCTGGCGGACGTGTCGGCCGTCA
TGGACAGCTTCTGACCATGATGGTGCCGGGGCGGCTACACGTCAAGCACCGCCTGGTGT
CTGATGTCAAGTGCCACCAAGATCCCGCACATCTGGCTCATGCTGAGCACCAAGATGCCTG
TCGTGTTTGACCGAAAGGCGTCGCGGGCTCACCAGGACTGGGCCCGGCTGCGCTGGTTTCG
TCACCATCCAGCCAGCCACATCGGAGCAGTATGAGTTGCGCTTCAGGCTGCTGGACCCGC
GGACACAGCAGGAGTGCGCCAGTGTGGCGTCATCCCCGTGGCTGCCTGCACCTTCGACG
TCCGAAACCTGCTGCCAACCGATCCTATAAGTTACCATCAAGAGGGCCGAGACCTCCA
CGCTGGTGTAAGAGCCCTGGAGGGACAGCCTCACCTGCACACCAAGCCGGAGCCCTGG
AGGGGCCCGCCCTCAGCCACTCTGTCTGAGAGATGATTTTCTAATATTTATCCACTAATA
AAGAAGAGTGTAATGCACATATGGAATTA
  
```


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for NM_024059 unedited <pre> NCCATGTTTCGANGATATTTGTNAATACGAACTTCACTATAGNNGGCGGCCGCAATTCC CGGGGTATCGTCGACCCACTGCGTCCGGGAGAGGCCCCAGCACTGACCTCCGTGGGGTG GAGATGAGGAGGATGAAAGGGTGTCTTCCTCCAGCATCTTCCTGAAGCTCCCGGATAAT GAGCACCACATGTGGCAGGCCTGGGCCTGGACAAGATGAAGCTGGGCAATCCCCAGTCCTT CCTGGACCAGGAGGAGGAGATGACCAGCAGCTGCTGGAACAGAGGCGTGGAAGACCTA CACCGAGCGCCGCAATGCCCTGCGTGAGTTCTGACCTCGGACCTGAGTCCGCACCTGCT CAAGCGCCACCACGCCCGCATGCAGCTGCTGCGTAAGTGCTCCTACTACATCGAGGTCCT GCCCAAGCACCTGGCCCTGGGCGACCAGAACCCGCTGGTGCTGCCTAGCGCCTTGTTCCT GCTCATCGACCCCTGGAAGTTCCAGCGCATGAAGAAGGTGGGCACAGCTCAGACCAAGAT CCAGCTCCTGCTGCTCGGGGACCTGTTGGAACAGCTCGACCATGGCCGTGCTGAGCTGGA TGCCCTGCTCCGGTCGCCAGACCACGGCCCTTCTGGCCGACTGGGCGCTGGTGGAGCG GCGGCTGGCGGACGTGTCGGCCGTCATGGACAGCTTCTGAACATGATGGTGCCGGGGCG GCTACACGTCAAGCACCCCTGGTGTCTGATGTCAGTGCCACCAAGATCCCGCACATCTG GCTCATGCTGAGCACCAAGATCCTGTCGTGTTGACCGAAAGCTTCGGCAGCTCACCAGA CTGGGCCCGGCTGCGCTGTTCTTACCATCCACCCGCCACATCAAGCCATTTGAGTTGCG CTTCAGCTG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_024059
Insert Size:	1124 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_024059.2</u> , <u>NP_076964.1</u>
RefSeq Size:	1124 bp
RefSeq ORF:	957 bp
Locus ID:	79025
UniProt ID:	<u>Q9BVV2</u>
Cytogenetics:	20q13.33