

Product datasheet for **SC123072**

TMEM60 (NM_032936) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM60 (NM_032936) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM60
Synonyms:	C7orf35; DC32
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF sequence for NM_032936 edited ATGAGAATGTCCTTGGCTCAGAGAGTACTACTCACCTGGCTTTTCACACTACTCTTCTTG ATCATGTTGGTGTGAACTGGATGAGAAAGCACCTTGGAACTGGTTCCATATTTATT CCAGTCTGGATATTTGATACTATCCTTCTTGTCTGCTGATTGTGAAAATGGCTGGGCGG TGTAAAGCTGGCTTTGACCCTCGACATGGATCACACAATATTAAGCAAGCTGGTA CCTCATTGCAATGTTACTTAAATTAGCCTTCTGCCTCGCACTCTGTGCTAACTGGAACA GTTTACTACCATGAATCTATCCTATGTCTTCATTCTTTATGGGCCTTGCTGGCTGGGC TTTAACAGAACTCGGATATAATGTCTTTTTGTGAGAGACTGA
5' Read Nucleotide Sequence:	>OriGene 5' read for NM_032936 unedited GACATCCTGGATTGTATCACACTCTCTATAGGCGGCCGACGCATTTCGACGGGCGGAAA AGNCAAGGATCCAACTGGCGAATTTGCTGATCTTCGCGTCCCTCTCCGCTTTCCGGCCG GCAGCGCTGCCAGGTATATTTCTTTTTTCCGATCCTGCAACAGCCTCTTTAACTGTT TAAATGAGAATGTCCTTGGCTCAGAGAGTACTACTCACCTGGCTTTTCACACTACTCTTC TTGATCATGTTGGTGTGAACTGGATGAGAAAGCACCTTGGAACTGGTTCCATATTT ATTCCAGTCTGGATATTTGATACTATCCTTCTTGTCTGCTGATTGTGAAAATGGCTGGG CGGTGTAAGTCTGGCTTTGACCCTCGACATGGATCACACAATATTAAGCAAGCTGGTA GTACCTCATTGCAATGTTACTTAAATTAGCCTTCTGCCTCGCACTCTGTGCTAACTGGA ACAGTTTACTACCATGAATCTATCCTATGTCTTCATTCTTTATGGGCCTTGCTGGCTGG GGCTTTAACAGAACTCGGATATAATGTCTTTTTGTGAGAGACTGACTTCTAAGTACATC ATCTCCTTTCTATTGCTGTTCAACAAGTTACCATTAAAGTGTCTGAATCTGTCAAGCTT CAGAATACCAGAGAACTGAGGAAAAATACCAAATGTAGTTTTATACTACTTCCATANAAC AGGATTGGNTGAATCACGGACTTCTAGTCAAACTACAGCTTAATTATTCAGCATTTGAGT TATTGAGATCCTNATTATCTCTTATGTAATAAAAGTTTGTGTTGGACCTCATTTTCTAC CTGAAAAAAAAAAAAAAAAAAAACTCGACTCTTGAT
Restriction Sites:	Please inquire


[View online »](#)

ACCN:	NM_032936
Insert Size:	809 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_032936.2</u> , <u>NP_116325.1</u>
RefSeq Size:	883 bp
RefSeq ORF:	402 bp
Locus ID:	85025
UniProt ID:	<u>Q9H2L4</u>
Cytogenetics:	7q11.23
Protein Families:	Transmembrane