

Product datasheet for **SC123488**

TRAPPC2B (BC032809) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRAPPC2B (BC032809) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRAPPC2B
Synonyms:	MIP-2A; SEDLP; SEDLP1; TRAPPC2P1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC032809 edited

```
GCAGAGGCGGTAGTGACACAGGCACAAGCTGACAGTGGCAGAAGCTCAGCTGACAAGGACT
GGGGACGGCGGTGCTTGTCTTGCCTTTGTGCGCCCCGCCCTCTCTCCCTGGCTGGA
CTTGCGGAGTCCCGGCCGAAGAAGCCGAGGAGCCATATATTGAAGACCATGTCTGGAAGC
TTCTACTTTGTAATTGTTGGTCACCATGATAATCCAGTTTTTGAAATGGAGTTTTTGCCA
GCTGGGAAGGCAGAATCCAAAGACGACCATCGTCATCTGAACCAGTTCATAGCTCATGCT
GCTCTCGACCTCGTAGATGAGAACATGTGGCTGTGCAACAACATGTACTTGAAAAGTGTG
GACAAGTTCAACGAGTGGTTTGTGTGAGCATTGTGACCGCGGGGCATATGAGGTTTATT
ATGCTTCATGACATAAGACAAGAAGATGGAATAAAGAAGTCTTTACTGATGTTTATGAT
TTATATATAAAATTTCAATGAATCCATTTTATGAACCAATTCTCTATTTCGATCAAGT
GCATTTGACAGAAAAGTTCAGTTTCTTGGGAAGAAACACCTTTTAAGCTGAATGGAGAAA
ATTCCAAAATAAATTATATCACCACAATGGTGATACTCAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAA
```



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC032809 unedited <pre> GTACGTCAGTATTTGTATACGACTCACTATAGGCGGCCGCGATTCCCGGGATGCAGAGGC GGTAGTGACACAGGCACAAGTACAGTGGCAGAAAGCTCAGCTGACAAGGACTGGGGACGG CGGTGTCCTTGTCTTGCCTTTGTGCGCCCCGCCCTCTCTCCCTGGCTGGACTTGCGGA GTCCCCGCCGAAGAACCCGAGGAGCCATATATTGAAGACCATGTCTGGAAGCTTCTACTT TGTAATTGTTGGTCACCATGATAATCCAGTTTTGAAATGGAGTTTTTCCAGCTGGGAA GGCAGAATCCAAAGACGACCATCGTCATCTGAACCAAGTTCATAGCTCATGCTGCTCTCGA CCTCGTAGATGAGAACATGTGGCTGTGCAACAACATGTACTTGAAAACGTGGACAAGTT CAACGAGTGGTTTGTGTGTCAGCATTGTGTCACCGCGGGGCATATGAGGTTTATTATGCTTCA TGACATAAGACAAGAAGATGGAATAAAGAACTTCTTTACTGATGTTTATGATTTATATAT AAAATTTTCAATGAATCCATTTTATGAACCAATTCTCTATTGATCAAGTGCATTTGA CAGAAAAGTTCAGTTTCTTGGGAAGAAACACCTTTTAAGCTGAATGGAGAAAATTCCAAA ATAAATTATATCACCCCAATGGGTGTATACTCAAAAAAGAGAGATCGGCCACGACATGG AAGTTAGAGTAGACCCAGGTCACTAAAAAGGCGGCCGCGGACATAACTGTTTCCTGG AAAAATCCCGGTGGGATTCTGTGGAACCTCACCAAGGGCTTCCCTGGCCCTTGGA GGTGCCACTTCGATTGCCACCAAGCCTGGCCGAAGAAGACTTAACGTGCTCCATTGTTN </pre>
Restriction Sites:	Please inquire
ACCN:	BC032809
Insert Size:	686 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC032809.1
RefSeq Size:	686 bp
Locus ID:	10597
Cytogenetics:	19q13.43