

Product datasheet for SC123500

LOC554207 (BC031469) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC554207 (BC031469) Human Untagged Clone
Tag:	Tag Free
Symbol:	LOC554207
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC031469 edited AAGAAGACTGCTTTTGGTGTTCCTTTAACCTGATGTCCCATTAATAATTTCAGAAGG AAAGGCCTGCTCATATAATTTCTACGGCCTCAGGCCAGAGACAGCCATACTTACAGGAC TGGAGGCTGCCACCACTGACATGTTCCACCAGATTCTTGTGGGCTCAAGAAGCATTCAA GCTTCATCCCCCTTCGTATTTATGAAATCCGGAGGTACTGGAGCAGCGCTGTATGTCAG CATCTGGCATTGTTCAATCAAGATGTTAGCTGGGACAGAAGACATAAGTCAGAACGCTGG AAGAACTGGGTCCCAATATTATAATAATCAATAAAGACAAAATATTTATAGTTTAA AAAAAAAAAA
5' Read Nucleotide Sequence:	>OriGene 5' read for BC031469 unedited NATTTTCCATTTGTCATACGACTCACTATAGGGCGGCCGCGGATCCCGGGGTATCGTCG ACCCACGCGTCCGAAGAAGACTGCTTTGGGTGTTTGCCTTAACCTGATGTCCCATTA ATAATTTCAGAAGGAAAGGCCTGCTCATATAATTTCTACGGCCTCAGGCCAGAGACAGCC ATACTTACAGGACTGGAGGCTGCCACCACTGACATGTTCCACCAGATTCTTGTGGGCTC AAGAAGCGATGCGAGCTTCATCCCCCTTCGTATTTATGAAATCCGGAGGTACTGGAGCAG CGCTGTATGTCCAGCATCTGGCATTGTTCAATCAAGATGTTAGCTGGGACAGAAGACATA AGTCAGAACGCTGGAAGAACTGGGTCCCAATATTATAATAATCAATAAAGACAAAATAT ATTTATAGTTTAAAAAAAAAAAAAAAAAGGGCGGCCGCGGTATAGCTGTTTCTGAACAGA TCCCGGGTGGCATCCCTGTGACCCCTCCCACTGCTCTCTGGCCCTGGAAGTTGCCAC TCCAGTGCCACCAAGCCTTGTCCTAATAAAATTAAGTTGCATCATTTTGTCTGACTAGGT GTCCTTCTATAATATTATGGGGTGGAGGGGGTGGTATGGGAACAGGGGGCAAGTTGGG AAACAACCTGGTAGGCCTGCGGGGGCTATTGGGAACCAACCTGGAGTGAATGGCACAAT CTTGGCTCACTGCAATCTCCGGCTCCTTGGTTTCAAGCGATTCTCTGGCCTCAGTCCC GAGTTGTTGGGGATTCCGGGCATGCGTGACCAGGCTCACCTAATTTTGTATTGGGA AGACGGGGTTTACCATATGGCCAGGCTGGTCTCCAATTCTATCTCAGGGATCTAACC ACCTTGGGCTTCAAATGCTGGAG
Restriction Sites:	Please inquire
ACCN:	BC031469


[View online »](#)

Insert Size:	371 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC031469.1</u> , <u>AAH31469.1</u>
RefSeq Size:	371 bp
Locus ID:	554207
Cytogenetics:	14q11.2