

Product datasheet for **SC123585**

RPL23AP82 (BC019348) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RPL23AP82 (BC019348) Human Untagged Clone
Tag:	Tag Free
Symbol:	RPL23AP82
Synonyms:	RPL23A_43_1761
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC019348 edited <pre> GGGTGTGGTACCGAGGCTTCAGCGGGTGCCGCCCGCCTAGAGGGAGTGGAGCGGTGAGCA CGTCAGGGGTGGGGGCGCAGGTCAAGCTTTCACCAAGTTTAAATCTTTGATGGGGTAA ATTTGAGCAATTTTCTCGACTTGTGACATTCGTTATTAAGTGAAGGAGGAGGAGAG GAACCCGGTCCTCTCCACACAGCCAGCAGAGAGCCTACGACTAGATTTGCATCTTTACG TCCTGCGCGGAGGCTGCTACACACATGCAGAAGTCATGCTGGTGGCCTGGACAGTGAAGG GAGAGAAGTGGATTTGGGAGACATTTAGGAGGAACAGTAAGAGGACCTTGTGCATGAATA ATTTGTTTCCACACTACAGAAAGAATCCCAGACTTCTTAGAGAACCCAGTGACTTTTGC ACCTTAAATCTGTGAAATCCTCATGTTTTCTTCTGCCGTATCCATAGTTCAAACAAAGAT GAGGCAAAGCTAGATGCATTCTGAAGGAACCAAGAAATTCCTCTTTATTTCTCTGG AATGAAATGAATTCTCTAGACCACAGTTCTAACCTTCAAAAACCAAACCTGTTTGTGAG ATCTCCTTCAAATACTACTGTAGACCCAGTGTTTATTCATTAATTTTAAATATTAA AAAAAAAAAAAAAAAAAAAA </pre>


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC019348 unedited <pre> NGGGGGTTCACATTTGTATACGACTCATATAGGGCGGCCGCAATTCGCACGAGGGGGTG TGGTACCGAGGCTTCAGCGGGTGCCGCCCGCTAGAGGGAGTGGAGCGGTGAGCACGTCA GGGGTGGGGGGCGCAGGTCAAGCTTTCACCAAGTTTAAATTCCTTGATGGGGTAAATTTG AGCAATTTTCTCGACTTGTGACATTCGTTATTAAGTACGAGGAATCAGGAGAGGAACC CGGTCCTCTCCACACAGCCAGCAGAGAGCCTACGACTAGATTTGCATCTTACGTCCTG CGCGGAGGCTGCTACACACATGCAGAAGTCATGCTGGTGGCCTGGACAGTGAAGGGAGAG AAGTGGATTTGGGAGACATTTAGGAGGAACAGTAAGAGGACCTTGTGCATGAATAATTTG TTTCCACACTACAGAAAGAATCCAGACTTCTTAGAGAACCCAGTGACTTTTGCACCTT ANATCTGTGAAATCCTCATGTTTTCTTCTGCCGTATCCATAGTTCANACAAAGATGAGGC AAAGCTAGATGCATTCTGAAGGAACCAAGAAATTCCTCTCTTATTTCTCTGGAATGA AATGAATTCTCTAGACCACCAAGTCTAACCTTCAAAAACCAACCTGTTTGTGAGATCTC CTTCAAATACTACTGTAGACCCAGTGTTTATTCATTAAATTTTAAAAATATTAATAA AAAAAAAAAAAACTCGACTCTAGATTGCCGCCGNGTCATAGCTGTTTCTGAACAGAT CCCGGGTGGCATCCCTGTGACCCCTCCCAAGTGCCTCTCTGGCCCTGGAAGTGCCACTC CAGTGCCCAACACCTTGTCTAATAAATTAAGTTTGCATCATTTGTCTGACTAGGTGTC CA </pre>
Restriction Sites:	Please inquire
ACCN:	BC019348
Insert Size:	678 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC019348.1
RefSeq Size:	678 bp
Locus ID:	284942
Cytogenetics:	22q13.33