

Product datasheet for **SC123592**

GAGE2E (BC018052) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GAGE2E (BC018052) Human Untagged Clone
Tag:	Tag Free
Symbol:	GAGE2E
Synonyms:	GAGE8; G antigen 2E
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC018052 edited GCTCACGTGAAGAACGCCAGGGAGCTGTGAGGCAGTGCTGTGGTTCTGCCGTCCGGA CTCTTTTCTCTACTGAGATTCATCTGTGTGAAATATGAGTTGGCGAGGAAGATCGACC TATCGGCCTAGACCAAGACGCTACGTAGAGCCTCCTGAAATGATTGGGCCTATGCGGCC GAGCAGTTCAGTGATGAAGTGGAACCAAGCAACCTGAAGAAGGGGAACCAAGCAACTCAA CGTCAGGATCCTGCAGCTGCTCAGGAGGGAGAGGATGAGGGAGCATCTGCAGGTCAAGGG CCGAAGCCTGAAGCTGATAGCCAGGAACAGGGTCACCCACAGACTGGGTGTGAGTGTGAA GATGGTCTGATGGGCAGGAGATGGACCCGCCAAATCCAGAGGAGGTGAAAACGCCTGAA GAAGGTGAAAAGCAATCACAGTGTTAAAAGAAGGCACGTTGAAATGATGCAGGCTGCTCC TATGTTGAAAATTGTTCAATAAATTCTCCAATAAAGCTTTACAGCCTTCTGCAAAAA AAAAAAAAAA


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC018052 unedited NACGGAGACTTTGATTTGTATACCATTACTATAGGCGGCCGCATAACTTCGTATAGCATA CATTACTACGAATTATGGCTCATGCCAAATCGGCCGAGCTCNGATTCTGTCGAGAGCGGGC TCACGCTGAANAACGCCAGGGAGCTGTGAGGCAGTGCTGTGTGGTTCTGCCGTCCGGAC TCTTTTTCCTCTACTGAGATTCTCTGTGTGAAATATGAGTTGGCGAGGAATATCGACCT ATCGGGGGAGAGCAAGACGCTACGTAGAGCCTCCTGAAATGATTGGGCCTATGCGGCCCG AGCAGTTCATTGATGAAGTGGAACCAACACCTGAGGAAGGGGAACCAGCAACTCAAC GTCAGGATCCTGCANCTGCTCATGAGGGAGAGGATGAGGGAGCATCTGCAGGTCAAGGGC CGAAGCCTGAAGCTGATAGCCAGGAACAGGGTCACCCACAGACTGGGTGTGAGTGTGAAG ATGGTCTCTGATGGGCAGGAGATGGACCCGCCAAATCCAGAGGAGGTGAAAACGCCTGAAG AAGGTGAAAAGCAATCACAGTGTTAAAAGAAGGCACGTTGAAATGATGCAGGCTGCTCCT ATGTTGGAAATTTGTTCAATAAATTCTCCAATAAAGCTTTACAGCCTTCTGCAAAAAA AAAAAAAAAACTCTCCAGCGCTGGATCCGGCCATAAGGGCCTGATCCTTCGAGGGGGGGC CCGGTACCGATATCAAGCTTGTGACTCTAGATTGCGGGCCGCGGGCCTAGCTGTTTCT GGACAGATCCCGGGGGGATCCCTGGGACCCCTCCCCAGTGCCTCTCTGCCCCGGGAAG TTGCCACTCCAGGGCCCACCAGCCTGGGCCTAATAAAATAAATTGGCTCG
Restriction Sites:	Please inquire
ACCN:	BC018052
Insert Size:	551 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC018052.1 , AAH18052.1
RefSeq Size:	551 bp
Locus ID:	26749
Cytogenetics:	Xp11.23