

Product datasheet for **SC123672**

UEVLD (BC011011) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UEVLD (BC011011) Human Untagged Clone
Tag:	Tag Free
Symbol:	UEVLD
Synonyms:	ATTP; FLJ11068; signaling molecule ATTP; ubiquitin-conjugating enzyme E2-like; UEV3; UEV and lactate/malate dehydrogenase domains
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC011011 edited

```

GGGGGACCTAGCTCGGGACCGCCTGGAGATGGAGTTCGACTGCGAGGGCCTGAGACGGC
TGCTTGCCAAGTACAAGTTCAGGGACCTAACTGTGGAAGAACTAAGGAATGTAATGTAT
TTTTCCACATTTCAAATATTCCATGGACACCTATGTTTTTAAAGATAGTTCTCAGAAAG
ACCTGCTGAATTTTACTGGCACAATTCCTGTGATGTATCAGGGTAATACATATAACATAC
CAATTCGTTTCTGGATTTTGGATTCTCACCTTTTCGCTCCCCCTATTTGCTTCTTGAAGC
CAACTGCAAAATATGGAATCTTAGTCGGAACATGTGGATGCTCAAGGCAGAATATATT
TGCCCTATCTCCAAACTGGAGCCATCCTAAATCTGTCATTGTTGGATTAATTAAGAAA
TGATTGCCAAGTTTCAAGAGGAACCTCCCATGTATTCTCTATCATCATCTGATGAGGCAC
GGCAGGTAGACTTGCTAGCCTATATTGCAAAAATCACTGAAGGTGTTTCAGATACAAATT
CAAAGAGCTGGGCAAATCATGAGAATAAAACAGTCAATAAAATTACTGTGGTTGGAGGTG
GAGAACTCGGTATTGCCTGCACATTAGCAATTTAGCAAAAGGTGGAGTCTCGCTCTGTTA
CCCAGGCTGGAGTGTAGTGGTACAATCTTGGCTCACTACAGCCTCCGCCTCCTGGGTTCA
AGCGATTCTCCTGCCTCAGCCTCCAGAGTAGCTGGGACTACAGGCATGCGCCACGACGCC
CAGCTAATTTTTGTATTTTAGCAGAGACAGGGTTTTGCCATGTTGGCCAGGCTGGTCTC
GAATCTTGACCTCAAGTGATCCACCTGCCCTGGCCTCCCAAAGTGCTGGGATTACAAGT
GTGAGCCACCGTACCTGGCCTGTATTACCTCTTAATTAGAAATTTGAGTTGTATTCAGTC
TTTTGTCACTAAATAAATGCTTCTATAAACATTCCAAAAA
AAAAA

```



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC011011 unedited <pre> GAAGTTCACATTTGTATACGACTCATATAGGCGGCACGCGAATTCGCCATTACGGCCGGG GGACCTAGCTCGGGACCGCCTAGGAGAGGAGTTCGACTGCGAGGGCCTGAGACGGCTGC TTGGCAAGTACAAGTTCAAGGACCTAACTGTGGAAGAACTAAAGAATGTAAATGTATTTT TCCCACATTTCAAATATTCATGGACACCTATGTTTTTAAAGATAGTTCTCAGAAAGACC TGCTGAATTTTACTGGCACAATTCCTGTGATGTATCAGGGTAATACATATAACATACCAA TTCGTTTTCTGGATTTTGGATTCTCACCTTTTCGCTCCCCCTATTTGCTTCTTGAAGCCAA CTGCAAATATGGGAATCTTAGTCGGAACATGTGGATGCTCAAGGCAGAAATATATTTGC CCTATCTCCAAAACCTGGAGCCATCCTAAATCTGTCAATTGTTGGATTAATTAAGAAATGA TTGCCAAGTTTCAAGAGGAACCTCCCATGTATTCTCTATCATCATCTGATGAGGCACGGC AGGTAGACTTGCTAGCCTATATTGCAAAAATCACTGAAGGTGTTTCAGATACAAATTCAG AGAGCTGGGCAAATCATGAGAATAAACAGTCAATAAAATTAAGTGGTTGGAGGTGGAG AACTCGGTATTGCCTGCACTTAGCAATTCAGCAAAAGTGGAGTCTCGCTCTGTTTACCC AGCTGGAGTGTAGTGGTACAATCTTGGCTCACTACAGCCTCCGCCTCTGGGTCAAGCG ATTCTCCTGCCTCAGCCTCCAGAATAGCTGGGACTACAGGCATGCGCCACGACGCCCAGC TAATTTTGTATTTTGTAGCAGAGACAGGGTTTGCCTCA </pre>
Restriction Sites:	Please inquire
ACCN:	BC011011
Insert Size:	1025 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC011011.1</u> , <u>AAH11011.2</u>
RefSeq Size:	1025 bp
Locus ID:	55293
Cytogenetics:	11p15.1
Gene Summary:	Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function]