

Product datasheet for **SC123713**

SNHG5 (BC009220) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SNHG5 (BC009220) Human Untagged Clone
Tag:	Tag Free
Symbol:	SNHG5
Synonyms:	bA33E24.2; C6orf160; LINC00044; NCRNA00044; U50HG
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC009220 edited <pre> GGCACGAGGGGGCGGTAGGTGAGTGGGTATTGCGGGCTAGTATCCGAGCAAAAGATGGTG GCGCAGGCCGAGTTAAGAGCTTTAATCCTGTGAAGACATCTTAGTGAAGAGTTTAGAGTG CTGAGAGTTGAAAAGCTTGACGTGGGAAACGTGCGGCCGGAAGTCCACATGTACTGAGGT TGAGTCGTGACGGCCACAGGCTCCGAGTTTTGGCGTGAGGAACCGCTGATCGGCCACGGG CGCCGAACTTGCTGGCCTCCGGCATGTGCCTGAGCGGCGGCGGAAAAACCACTTAATTG GGGCGGAGGGTTAGTTTTAACAGCAAAGGGCCTTTACTAAAAATGGCGAAGGCCTTCCGTC GGCGTTGTTTTAAATGGGAAGCCTCGACCCTGTATTGAACTGAGCTGTTTGAAGGCGG CGTTGTGTGCAATTCGATTAATGAAGGGGAAGGGTTTTGTGTGAAAAACGCCTTGGAG TGTGACATTTCTGCGAGAATGCTTAAATACCGATTTCCCGCAGGAACAATGGCGCTGTCT TCAGTGGCACAGTGGAGCAGCTCTGAAGATGCAAAGATACACGAAAAAATTCCAGAACA TCTGGGAGAATATTTAATGGAAAATCGCTTGTTAAACCTGACACTTTTAACAGTGAAC AGCGTTCTGAGTGTGGACGAGTAGCCAGTGAAGATAATGAATGTGCAATGTGACTGACTA GCAGCTTCATTTTGAATGAGGTCGCTGTCTGCCATTGATAGAGGCCAGATTGTCTTGG AAGTTCCAAAGTTGCAACGATTTCTGGCTAGTGCCACGAGGTTTACTTGACTGTTGTGTG AAAAGCTGATAAGAAAACCATCCAGAAAAAAGCTCTTCGTTTTACAAACATGAAAAATAA ACATGTAATTTTGGATTA A </pre>


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC009220 unedited <pre> NNNNTTAGTCAAAATTGTATACGACTCACTATAGGCGGCCCGCAATTCGCACGAGGGGC GGTAGGTGAGTGGGTATTGCGGGCTAGTATCCGAGCAAACTGGTGGCGCAGGCCGAGTT AAGAGCTTTAATCCTGTGAAGACATCTTAGTGAAGAGTTTAGAGTGCTGAGAGTTGAAAG CTTGACAGTGGGAAACGTGCGGCCGACTGCCACATGTACTGAGGTTGAGTCGTGACGGC CACAGGCTCCGAGTTTGGCGTGAGGAACCGCTGATCGGCCACGGGCGCCGAACCTTGCTG GCCTCCGCGCATGTGCCTGAGCGGCGCGGAAAAACACCTTAATTGGGGCGGAGGGTTAG TTTTAACAGCAAAGGGCCTTTACTAAAAATGGCGAATGCCTCCGTCGGCGTTGTTTTAAA ATGGGAAGCCTCGACCCTGTATTGAAACTGAGCTGTTTGAAGGCGGCGTTGTGTGAATT CGGATTAATGAAGGGGAAGGTTTTGTGTGAAAAACGCCTTGGAGTGTGACATTCTGC GAGAATGCTTAAATACCGATTTCCCGCAGGAACAATGGCGCTGTCTTCAGTGGCAGAGT GAGCAGCTCTGAAGATGCAAAGATACACGAAAAAACTTCCAGAACATCTGGGAGAATATT TAATGGAAATCGCTTGGTTAAACCTGACACTTTTAACAGTGAACAGCGTTCTGAGTGT GGACGAGTAGCCAGTGAAGATAATGAATGNNTCGATGTGACTGACTAGCAGCTTCATTTT GAATGAGGGTCGCTGTCTGCCCATGATAGAGCCAGATTGTCTTTGGAGTTCAAAGTTGC ACGATTTCTGGCTAGTGCCACGAGGTTACTTGACTGNTGTGTGAAAGCTGATAGAAAACC ATCCAGAAA </pre>
Restriction Sites:	Please inquire
ACCN:	BC009220
Insert Size:	961 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC009220.1 , AAH09220.1
RefSeq Size:	961 bp
Locus ID:	441164
Cytogenetics:	6q14.3
Gene Summary:	DISCONTINUED: This record has been withdrawn by NCBI because the model on which it was based was not predicted in a later annotation.